


Food, Feed, Fermentation, Biofuels, Wine, Brewing & Dairy Industries





# Megazyme Test Kits and Reagents

## **Purity. Quality. Innovation.**

**Barry V. McCleary, PhD, DScAgr**

*Innovative test methods with exceptional technical support and customer service.*  
**The Megazyme Promise.**

Megazyme was founded in 1988 with the specific aim of developing and supplying innovative test kits and reagents for the cereals, food, feed and fermentation industries. There is a clear need for good, validated methods for the measurement of the polysaccharides and enzymes that affect the quality of plant products from the farm gate to the final food.

The commitment of Megazyme to “Setting New Standards in Test Technology” has been continually recognised over the years, with Megazyme and founder, Dr McCleary receiving a number of business and scientific awards.

Today, Megazyme is a recognised world leader in the development of high quality, innovative test technology for many industries and offers the complete toolbox for all aspects of carbohydrate research. Over 500 products are available, 80% of which are unique to Megazyme and these are spread across five main categories; Assay kits, Enzymes, Enzyme Substrates, Carbohydrates and Equipment.



Megazyme has its research and manufacturing facilities in Bray, Ireland. Exports account for over 98% of sales with Megazyme currently shipping to 90 countries spread across every continent. To provide the best customer experience possible, we offer worldwide express delivery with time-stamped tracking on all of our orders and any technical queries on our products are answered by our scientist within 48 hours.

All products can be ordered at [www.megazyme.com](http://www.megazyme.com) and our website also acts as the ultimate information resource for all Megazyme products. Customers can download a data booklet for every assay kit that outlines clear step-by-step instructions for their use. Data booklets for other product classes contain detailed physico-chemical information from pH and temperature stability ranges for enzymes to NMR spectra and HPLC/IC chromatograms of enzyme substrates or carbohydrates. Product specific Mega-Calc™ Excel™ based calculation tools, certificates of analysis (COAs), material safety data sheets (MSDS) and frequently asked questions (FAQs) can also be found at [www.megazyme.com](http://www.megazyme.com).

**Professor Barry V. McCleary**  
Chief Executive Officer

# Table of Contents



## Assay Kits

Introduction	5
All Assay Kits	9
Industry Specific Assay Kits	10
Individual Assay Kit Descriptions	17

## Enzyme Substrates

Insoluble Chromogenic Substrates	67
Soluble Chromogenic Substrates	68
Enzyme Tablet Tests	70
Colourimetric Oligosaccharides	71

## Enzymes

64

## Carbohydrates

Oligosaccharides	74
Polysaccharides	77

## Equipment

79

Kit	Cat. No.	Page	Kit	Cat. No.	Page
Acetaldehyde	K-ACHYD	17	L-Glutamine / Ammonia (Rapid)	K-GLNAM	39
Acetic Acid (ACS; Analyser format)	K-ACETAF	17	Glycerol	K-GCROL	40
Acetic Acid (ACS; Manual format)	K-ACET	18	Glycerol (ADP-GK format)	K-GCROLGK	40
Acetic Acid (AK; Analyser format)	K-ACETAK	18	D-3-Hydroxybutyric Acid	K-HDBA	41
Acetic Acid (ADP-GK format)	K-ACETGK	19	myo-Inositol	K-INOSL	41
Acetic Acid (RM; Rapid manual format)	K-ACETRM	19	D-Isocitric Acid	K-ISOC	42
Ammonia (Rapid)	K-AMIAR	20	D-Lactic Acid	K-DATE	42
$\alpha$ -Amylase ("Ceralpha" method)	K-CERA	20	D-/L-Lactic Acid	K-DLATE	43
$\alpha$ -Amylase ("Sprout damaged" method)	K-AMYLSL	21	L-Lactic Acid	K-LATE	43
$\beta$ -Amylase ("Betamyl-3" method)	K-BETA3	21	Lactose / D-Galactose (Rapid)	K-LACGAR	44
Amylose / Amylopectin	K-AMYL	22	Lactose / Sucrose / D-Glucose	K-LACSU	44
Amyloglucosidase / Glucoamylase	K-AMG	22	Lactulose	K-LACTUL	45
Arabinan	K-ARAB	23	D-Malic Acid	K-DMAL	45
L-Arabinose / D-Galactose (Rapid)	K-ARGA	23	L-Malic Acid	K-LMAL	46
L-Arginine / Urea / Ammonia (Rapid)	K-LARGE	24	L-Malic Acid (Analyser format)	K-LMALAF	46
L-Ascorbic Acid	K-ASCO	24	L-Malic Acid (Liquid Ready Reagents)	K-LMALQR	47
L-Asparagine / L-Glutamine / Ammonia (Rapid)	K-ASNAM	25	L-Malic Acid (MegaQuant™ format)	K-LMALMQ	47
Aspartame	K-ASPTM	25	Malt Amylase	K-MALTA	48
Available Carbohydrates / Dietary Fiber	K-ACHDF	26	Maltose / Sucrose / D-Glucose	K-MASUG	48
endo-1,4- $\beta$ -Glucanase (cellulase)	K-CELLG3	26	D-Mannitol / L-Arabitol	K-MANOL	49
Citric Acid	K-CITR	27	D-Mannose / D-Fructose / D-Glucose	K-MANGL	49
Ethanol	K-ETOH	27	Pectin Identification	K-PECID	50
Fiber (Total Dietary Fiber)	K-TDFR	28	Phytic Acid (Total Phosphorus)	K-PHYT	50
Fiber (Integrated Total Dietary Fiber)	K-INTDF	28	Primary Amino Nitrogen (NOPA)	K-PANOPA	51
Formic Acid	K-FORM	29	Pyruvic Acid	K-PYRUJ	51
Fructan (Hexokinase format)	K-FRUCHK	29	Raffinose / D-Galactose	K-RAFGA	52
Fructan (PAHBAH format)	K-FRUC	30	Raffinose / Sucrose / D-Glucose	K-RAFGA	52
D-Fructose / D-Glucose	K-FRUGL	30	L-Rhamnose	K-RHAMNOSE	53
D-Fructose / D-Glucose (Liquid Ready Reagents)	K-FRGLQR	31	D-Sorbitol / Xylitol	K-SORB	53
D-Fructose / D-Glucose (MegaQuant™ format)	K-FRGLMQ	31	Starch Damage	K-SDAM	54
L-Fucose	K-FUCOSE	32	Starch (Resistant Starch)	K-RSTAR	54
Galactomannan	K-GALM	32	Starch (Total Starch; GOPOD format)	K-TSTA	55
Beta-Glucan (Mixed linkage)	K-BGLU	33	Starch (Total Starch; Hexokinase format)	K-TSHK	55
Beta-Glucan (Yeast and mushroom)	K-YBGL	33	Succinic Acid	K-SUCC	56
Beta-Glucan (Yeast-Enzymatic)	K-EBHLG	34	Sucrose / D-Fructose / D-Glucose	K-SUFRG	56
Beta-Glucanase (Malt and microbial)	K-MBGL	34	Sucrose / D-Glucose (GOPOD format)	K-SUCGL	57
Glucomannan	K-GLUM	35	Sulphite (Total SO <sub>2</sub> )	K-TSULPH	57
D-Gluconic Acid / D-Glucono- $\delta$ -lactone	K-GATE	35	Sulphite (Total SO <sub>2</sub> ; Enzymatic)	K-ETSULPH	58
D-Glucosamine	K-GAMINE	36	Sulphite (Total & Free SO <sub>2</sub> )	K-SULPH	58
D-Glucose (GOPOD format)	K-GLUC	36	Tartaric Acid	K-TART	59
D-Glucose (Hexokinase format)	K-GLUHK	37	Trehalose	K-TREH	59
Glucose Oxidase	K-GLOX	37	Urea / Ammonia (Rapid)	K-URAMR	60
D-Glucuronic Acid / D-Galacturonic Acid	K-URONIC	38	Xylanase (Azo-Wax format)	K-AZOWAX	60
$\alpha$ -Glucuronidase	K-AGLUA	38	Xylanase (Xylazyme AX format)	K-XYLS	61
L-Glutamic Acid (MSG)	K-GLUT	39	D-Xylose	K-XYLOSE	61
			Limit-Dextrinase / Pullulanase	K-LDPU	62



# ASSAY KITS...





## Principles of Test Procedures

In general terms, enzymes are catalytic proteins that convert one compound into another, and such reactions frequently occur without any visible sign that they have taken place. However, certain enzymatic reactions result in either a “colour” being produced or consumed, and the intensity of the colour change can be measured using a common spectrophotometer.

Invented in the 1950s, the spectrophotometer in its various guises is today one of the most commonly used analytical instruments. The spectrophotometer is a powerful analytical instrument because it can measure changes in absorbance very accurately and quickly. The enzymatic analysis reaction itself is performed in a plastic or glass cuvette that sits between the source of light and the detector inside the spectrophotometer (as depicted simplistically in figure 1). A known amount of light travels through the cuvette and the amount that emerges is quantified by the detector. The change in intensity as the light passes through the reaction solution in the cuvette is recorded as an absorbance reading.




Figure 1. Simplistic representation of enzymatic analysis being conducted using a spectrophotometer.

Modern enzymatic analysis test kits generally contain all reagents necessary to perform the assay, e.g. buffer, cofactor(s), trigger enzyme and standard solution, in an easy to use form that is stable for > 2 years (even while in use). During a typical enzymatic analysis, deionised water is mixed with buffer, cofactor and sample. Then an absorbance reading ( $A_1$ ) is taken just before addition of the “trigger” enzyme (specific for the analyte in question), after which the reaction takes place (see figure 2). When the reaction has finished (i.e. the “endpoint” has been achieved), a second absorbance reading is taken ( $A_2$ ). The difference between these two absorbance readings (i.e.  $A_1 - A_2$ ) is called the change in absorbance (or  $\Delta A$ ) and is directly related to analyte content. It is this value, after correction with a blank reading (reaction containing no sample), that is used to calculate the concentration of the analyte in the sample (typically as g/L), using a simple factor (e.g. 0.2535 in the case of the the acetic acid AK / PTA format). A sample calculation is shown below.

### Typical acetic acid calculation

$$\begin{aligned}
 A_1 (\text{blank}) &= 1.400 & A_2 (\text{blank}) &= 1.398 \\
 A_1 (\text{sample}) &= 1.420 & A_2 (\text{sample}) &= 0.650 \\
 \Delta A_{\text{acetic acid}} &= (A_1 (\text{sample}) - A_2 (\text{sample})) - (A_1 (\text{blank}) - A_2 (\text{blank})) \\
 \Delta A_{\text{acetic acid}} &= (1.420 - 0.650) - (1.400 - 1.398) \\
 &= 0.768 \\
 \text{Thus the concentration of acetic acid} &= 0.768 \times 0.2535 \\
 &= 0.1947 \text{ g/L}
 \end{aligned}$$

When using Megazyme test kits, calculations can either be performed manually as illustrated above or by using a free Excel™-based calculator, called “Mega-Calc™”, downloadable from the product page on the Megazyme website ([www.megazyme.com](http://www.megazyme.com)).




Figure 2. Decrease in absorbance at 340 nm on incubation of 0-35 µg of acetic acid with acetate kinase in the acetic acid AK / PTA format.

## Types of Enzymatic Assays

### Enzymatic Assays Based on $\text{NAD}^+$ / $\text{NADH}$ / $\text{NADP}^+$ / $\text{NADPH}$

Megazyme assay kits usually involve enzymes that either directly, or indirectly (via other enzymes), produce or consume a compound called NADH (or NADPH), that although invisible to the human eye, absorbs light strongly at a wavelength of 340 nm (extinction coefficient  $[\epsilon] = 6300 \text{ M}^{-1} \text{ cm}^{-1}$ ). Figure 3 depicts the various types of enzymatic reactions, either consuming or producing NADH (or NADPH), that are employed in many enzymatic test kits (Reaction 1). As can be seen from Figure 3, it is sometimes necessary to include an additional reaction in order to obtain quantitative results (Reaction 2).

This reaction is catalysed by the enzyme diaphorase, in the presence of a compound called INT, which converts the NADH (or NADPH) produced in the first reaction into a red coloured compound called INT-formazan. The resulting change in absorbance at 492 nm can be used to quantify the analyte of interest as discussed previously.




Figure 3. Enzymatic Assays Based on  $\text{NAD}^+$  /  $\text{NADH}$  /  $\text{NADP}^+$  /  $\text{NADPH}$

### Enzymatic Assays Based on GOPOD determination

A number of Megazyme test procedures are based on the ability to quantify glucose using the Megazyme Glucose oxidase/Peroxidase system (GOPOD). This can be applied directly to measure glucose in a sample (K-GLUC) but also finds use in the measurement of analytes that can be stoichiometrically converted to glucose (e.g. K-LACSU). The principle of the GOPOD system is shown in Figure 4. In Reaction 1, glucose is converted to glucono- $\delta$ -lactone by glucose oxidase with the production of  $\text{H}_2\text{O}_2$ . In Reaction 2,  $\text{H}_2\text{O}_2$  is used by peroxidase to form a quinoneimine that absorbs at 510 nm. The resulting change in absorbance at 510 nm can be used to quantify the analyte of interest as discussed previously.




Figure 4. Enzymatic Assays Based on GOPOD determination

### Enzymatic Assays Based on Chromogenic Substrates / Colourimetric Oligosaccharides

Certain assay kits involving the measurement of hydrolytic enzyme activity are based on the use of chromogenic substrates (e.g. K-AZOWAX) or colourimetric oligosaccharides (e.g. K-CELLG3). Chromogenic substrates are dyed and cross-linked polysaccharides that release colour upon hydrolysis. The colour released is directly proportional to the enzymatic activity being analysed. Colourimetric oligosaccharides are defined oligosaccharides covalently bound to chromogenic or fluorogenic moiety. These moieties are cleaved by hydrolytic enzymes releasing the desired chromophore or fluorophore which can be measured using a spectrophotometer. The principles of these assays are discussed in detail on page 71.

## Complete Solutions for Biochemical Analysis

Megazyme has teamed up with Awareness Technology to offer complete analytical solutions in the form of Megazyme assay kits (assay reagents) and associated equipment required to perform the analyses. Suitable assay formats and equipment combinations are available for customers at every level, from high throughput, large analytical laboratories right down to the smallest analytical facilities having minimal scientific equipment.



### Assay kits

While many of the Megazyme assay kits for bio-analysis are based on Boehringer Mannheim Methods of Analysis, continual innovation has enabled the introduction of new developments and improvements. Advanced enzymes offer reduced reaction times, and/or increased stability or alternative biochemical reactions that have been evaluated and implemented. Cofactor stability is enhanced by offering these in a stable tablet form or as improved formulations that provide extended stability in solution. Traditionally, assay kits for bio-analysis were supplied for use in the manual spectrophotometer format. Megazyme has adapted this to offer either a specific assay kit or modified procedure providing a comprehensive range of assays that can be used across the various assay formats:

- **Manual format** - manual format for use with cuvettes and a standard UV/Vis spectrophotometer.
- **Microplate format** - automated or manual format for use with standard 96-well microplates and microplate reader.
- **Auto-analyser format** - automated formats for use with any auto-analyser model.

Examples of the available assay kits covering the various assay formats are outlined below for the analysis of **L-malic acid**, a key parameter measured in wine analysis:

1. **K-LMALAF** is an advanced formulation of the traditional kits offered by other manufacturers for use in the auto-analyser assay format. The advantages of this kit include (i) PVP is incorporated to prevent sample interference (e.g. from tannins in wine), (ii) very stable R1 and R2 reagents, (iii) all reagents are stable for > 2 years during use (both enzymes are supplied as ammonium sulphate suspensions), and (iv) linear calibration ( $R^2 \sim 0.9994$ ) up to 80  $\mu\text{g/mL}$  in final reaction solution (see page 46).
2. **K-LMALQR** is supplied "ready to use" as liquid stable formulations recommended for high throughput use in the microplate or auto-analyser assay formats. The advantages of this kit include (i) no reagent preparation required (ii) very simple format, (iii) very cost effective (iv) PVP is incorporated to prevent sample interference (e.g. from tannins in wine) (iv) reagents are stable for > 2 years (see page 47).
3. **K-LMAL** is an advanced formulation of the traditional kits offered by other manufacturers for use in the manual assay format. The advantages of this kit include (i) PVP is incorporated to prevent sample interference (e.g. from tannins in wine), (ii) all reagents are stable for > 2 years during use (both enzymes are supplied as ammonium sulphate suspensions), (iii) very rapid reaction ( $\sim 3$  min) and (iv) an Excel™ based calculator (Mega-Calc™) is available on-line for hassle-free raw data processing (see page 46).
4. **K-LMALMQ** is a novel **manual assay format** used with the **MegaQuant™ Colorimeter** and is recommended for users who do not possess a laboratory and / or analytical expertise. The advantages of this kit include (i) samples do not need to be sent out for contract analysis, (ii) no spectrophotometer is required, (iii) very simple format, (iv) very cost effective, and (v) accurate and reliable with all samples (including red wine) (see page 47).



## Equipment

- The **ChemWell® 2910** is a fully automated open system analyser with optimised programmed protocols ready for use with Megazyme assay kits in **auto-analyser assay format**. It is capable of performing 200 biochemistry assays per hour in micro-well strips using a standard 96-well layout. On completion of assays the micro-wells can be automatically washed for re-use and, being equipped with 27 reagent positions and 96 sample positions, this provides the capability to perform continuous tests without intervention.
- The **ChemWell®-T** is a fully automated open system analyser with optimised programmed protocols ready for use with Megazyme assay kits in **auto-analyser assay format**. It is capable of performing 100 biochemistry assays per hour in 1 cm pathlength cuvette strips. Being equipped with 35 combined reagent and sample positions, it is capable of performing 40 tests without intervention. Cuvettes are easy to load and unload and the “continuous loading” capability allows tests to continue with ease.
- Ideally suited to lower throughput users, the **Stat Fax® 4500 Chemistry Analyser** is a compact, economical, standalone spectrophotometer supplied with programmed protocols for an extensive range of the Megazyme assay kits in **manual assay format**. The protocols offer ease of use for the analyst, providing on-board step-by-step instructions and automated calculation of results. In addition, the **Stat Fax® 4500** can be used as a standard spectrophotometer in various modes offering even greater flexibility.
- The **MegaQuant™ Colorimeter** is an extremely affordable alternative that can be used to measure residual sugar (D-fructose and D-glucose) and L-malic acid in wine. This is a basic, handheld colorimeter that can provide access to the measurement of these key parameters of wine quality for those without any laboratory equipment.

## Method Validation and Recognition by Official Bodies

Many of the analytical methods developed by Megazyme have become the official reference method recommended by the appropriate regulatory body for the area in question. In the area of dietary fiber for example, the international body that sets guidelines for national governments is called CODEX Alimentarius. Of the 14 methods currently approved by CODEX for the measurement of total dietary fiber and dietary fiber components, four were developed by Megazyme.

The table below shows the official bodies that have approved various Megazyme test procedures. Also shown are their abbreviations which are used from pages 16-62 in the relevant individual assay descriptions.

<b>AIJN</b>	Association of the Industry of Juices and Nectars from Fruits and Vegetables	<b>IFU</b>	International Federation of Fruit Juice Producers
<b>AACC</b>	American Association of Cereal Chemists	<b>IOCCC</b>	Office International du Cacao, du Chocolat et de la Confiserie
<b>AOAC</b>	Association of Official Analytical Chemists	<b>ISO</b>	International Standard Organisation
<b>ASBC</b>	American Society of Brewing Chemists	<b>IUPAK</b>	International Union of Pure and Applied Chemistry
<b>CCFRA</b>	Campden & Chorleywood Food Research Association Group	<b>JECFA</b>	Joint FAO / WHO Expert Committee on Food Additives
<b>CODEX</b>	CODEX Alimentarius	<b>MEBAK</b>	Central European Brewing Committee for Analysis
<b>DIN</b>	Deutsche Industrie Norm (German Standard)	<b>NBN</b>	Norme Belge (Belgian Standard)
<b>EBC</b>	European Brewery Convention	<b>NEN</b>	Nederlandske Norm (Dutch Standard)
<b>EEC</b>	Council of European Communities	<b>NF</b>	Normes Françaises (French Standard)
<b>EN</b>	European Norms	<b>NMKL</b>	Nordisk Metodikkomité for Næringsmidler (Nordic Committee of Food Analysis)
<b>GOST</b>	GOSSTANDART (State Committee of the Russian Federation for Standardisation and Metrology)	<b>OIV</b>	Office International de la Vigne et du Vin (International Wine Office)
<b>ICC</b>	International Association for Cereal Science and Technology	<b>RACI</b>	Royal Australian Chemical Institute
<b>ICUMSA</b>	International Commission for Uniform Methods of Sugar Analysis	<b>UKMBI</b>	United Kingdom Milling and Baking Industries
<b>IDF</b>	International Dairy Federation		

## ASSAY KITS

Cat. No.	Product
K-ACHYD	Acetaldehyde
K-ACETAF	Acetic Acid (ACS; analyser format)
K-ACET	Acetic Acid (ACS; manual format)
K-ACETAK	Acetic Acid (AK; analyser format)
K-ACETGK	Acetic Acid (ADP-GK format) <sup>NEW</sup>
K-ACETRM	Acetic Acid (AK; rapid manual format)
K-AMIAR	Ammonia (Rapid)
K-CERA	$\alpha$ -Amylase (Ceralpha method)
K-AMYLSD	$\alpha$ -Amylase (Sprout damaged) <sup>NEW</sup>
R-AMHR4	Amylase HR Reagent - 4 Vials
R-CAAR4	Ceralpha; $\alpha$ -Amylase Reagent - 4 Vials
K-BETA3	$\beta$ -Amylase (Betamyl-3 method)
R-BAMR3	Betamyl-3; $\beta$ -Amylase Assay Reagent - 4 Vials
K-AMYL	Amylose/Amylopectin
K-AMG	Amyloglucosidase <sup>NEW</sup>
R-AMGR3	Amyloglucosidase Assay Reagent - 4 Vials
K-ARAB	Arabinan
K-ARGA	L-Arabinose / D-Galactose (Rapid) <sup>NEW</sup>
K-LARGE	L-Arginine / Urea / Ammonia (Rapid)
K-ASCO	L-Ascorbic Acid (L-Ascorbate)
K-ASNAM	L-Asparagine / L-Glutamine / Ammonia (Rapid)
K-ASPTM	Aspartame
K-ACHDF	Available Carbohydrates / Dietary Fiber
K-CELLG3	endo-Cellulase <sup>NEW</sup>
R-CELLFLR	CellafLOUR (Cellulase Assay Reagent) <sup>NEW</sup>
K-CITR	Citric Acid
K-ETOH	Ethanol
K-TDFR	Fiber (Total Dietary Fiber)
K-INTDF	Fiber (Integrated Total Dietary Fiber) <sup>NEW</sup>
K-TDFC	Fiber Controls (Total Dietary Fiber)
K-FORM	Formic Acid
K-FRUCHK	Fructan (Hexokinase format)
K-FRUC	Fructan (PAHBAH format)
K-FRUGL	D-Fructose / D-Glucose
K-FRGLQR	D-Fructose / D-Glucose (Liquid Ready Reagents)
K-FRGLMQ	D-Fructose / D-Glucose (MegaQuant™ format)
K-FUCOSE	L-Fucose
K-GALM	Galactomannan (Carob or guar)
K-BGLU	$\beta$ -Glucan (Barley & oat; mixed linkage)
K-YBGL	$\beta$ -Glucan (Yeast & mushroom)
K-EBHLG	$\beta$ -Glucan (Yeast-enzymatic)
K-MBGL	$\beta$ -Glucanase (Malt & microbial)
K-GLUM	Glucomannan
K-GATE	D-Gluconate / D-Glucono- $\delta$ -lactone
K-GAMINE	D-Glucosamine (D-Glucosamine sulphate)
K-GLUC	D-Glucose (GOPOD format)
R-GLC4	Glucose Determination Reagent
K-GLUHK	D-Glucose (Hexokinase format)
K-GLOX	Glucose Oxidase
K-URONIC	D-Glucuronic Acid/D-Galacturonic Acid

Cat. No.	Product
K-AGLUA	$\alpha$ -Glucuronidase <sup>NEW</sup>
K-GLUT	L-Glutamic Acid (MSG)
K-GLNAM	L-Glutamine / Ammonia (Rapid)
K-GCROL	Glycerol
K-GCROLGK	Glycerol (ADP-GK format)
K-HDBA	D-3-Hydroxybutyric Acid
K-INOSL	myo-Inositol <sup>NEW</sup>
K-ISOC	D-Isocitric Acid
K-DATE	D-Lactic Acid
K-DLATE	D- / L-Lactic Acid
K-LATE	L-Lactic Acid
K-LACGAR	Lactose / D-Galactose (Rapid)
K-LACSU	Lactose / Sucrose / D-Glucose
K-LACTUL	Lactulose
K-LDPU	Limit-Dextrinase / Pullulanase <sup>AVAILABLE SOON</sup>
K-DMAL	D-Malic Acid
K-LMAL	L-Malic Acid
K-LMALAF	L-Malic Acid (Analyser format)
K-LMALQR	L-Malic Acid (Liquid Ready Reagents)
K-LMALMQ	L-Malic Acid (MegaQuant™ format)
K-MALTA	Malt Amylase
K-MASUG	Maltose / Sucrose / D-Glucose
K-MANOL	D-Mannitol / L-Arabitol
K-MANGL	D-Mannose / D-Fructose / D-Glucose
K-PECID	Pectin Identification
K-PHYT	Phytic Acid (Total Phosphorus)
K-PANOPA	Primary Amino Nitrogen (NOPA)
K-PYRUV	Pyruvic Acid
K-RAFGA	Raffinose / D-Galactose
K-RAFGL	Raffinose / Sucrose / D-Glucose
K-RHAMNOSE	L-Rhamnose <sup>NEW</sup>
K-SORB	D-Sorbitol / Xylitol
K-SDAM	Starch Damage
K-RSTAR	Starch (Resistant Starch)
K-RSTCL	Starch Controls (Resistant Starch)
K-TSTA	Starch (Total Starch; GOPOD format)
K-TSHK	Starch (Total Starch; Hexokinase format) <sup>NEW</sup>
K-TSCK	Starch Controls (Total Starch)
K-SUCC	Succinic Acid
K-SUFRG	Sucrose / D-Fructose / D-Glucose
K-SUCGL	Sucrose / D-Glucose
K-TSULPH	Sulphite (Total SO <sub>2</sub> )
K-ETSULP	Sulphite (Total SO <sub>2</sub> ; Enzymatic)
K-SULPH	Sulphite (Total & Free SO <sub>2</sub> )
K-TART	Tartaric Acid
K-TREH	Trehalose
K-URAMR	Urea / Ammonia (Rapid)
K-AZOWAX	Xylanase (Azo-Wax format)
K-XYLS	Xylanase (Xylazyme AX format)
K-XYLOSE	D-Xylose



## Food Industry Test Kits

Analyte	Cat. No.	Analyte Significance	Advantages of Megazyme Test Kits
<b>Acetic Acid</b>	K-ACETRM K-ACETAK K-ACETGK	Common food component	K-ACETAK (auto) / K-ACETRM (manual) are very rapid acetate kinase (AK) based kits with excellent linearity. K-ACETGK is a new rapid, auto-analyser assay kit employing AK and phosphotransacetylase. Stable reagents
<b>Ammonia</b>	K-AMIAR	Common food component	K-AMIAR has a very rapid reaction rate (~ 3 min at room temperature). Manual and auto-analyser applications
<b>Amylose / Amylopectin</b>	K-AMYL	Ratio of these components affects the rate of digestion and utilisation of starch	Novel kit, stable reagents
<b>L-Asparagine / L-Glutamine / Ammonia</b>	K-ASNAM	Acrylamide precursors in the production of fried, roasted, toasted potato or other food products	Novel product, enabling all three analytes to be determined in less than 20 min. Manual and microplate format procedures given
<b>L-Ascorbic Acid</b>	K-ASCO	Naturally found in fruits and vegetables, or supplemented in processed foods	Rapid reaction, stable reagents
<b>Available Carbohydrates / Dietary Fiber</b>	K-ACHDF	Sugars rapidly digested and absorbed, and dietary fibre	Novel procedure, stable reagents
<b>β-Glucan (Mixed linkage)</b>	K-BGLU	Major cell-wall polysaccharide of barley and oats	Rapid reaction, stable reagents, only enzymatic kit available. AOAC Method 995.16; AACC Method 32-23.01; ICC Standard No. 166; RACI Standard Method
<b>Citric Acid</b>	K-CITR	Common food component / additive	Ideal for manual and auto-analyser applications
<b>Ethanol</b>	K-ETOH	Found in small amounts in many foods	Rapid reaction, stable reagents (AIDH supplied as a stable suspension)
<b>Fructan</b>	K-FRUC K-FRUCHK	Common component in many foods such as onions and seeds	Novel assays, rapid reaction, stable reagents; AOAC Method 999.03; AACC Method 32-32.01
<b>D-Fructose / D-Glucose</b>	K-FRUGL K-FRGLMQ K-FRGLQR	Very common food sugars, e.g. from high fructose corn syrup supplementation	Ideal for manual and auto-analyser applications. Stable reagents. Choice of spectrophotometric or simple colorimeter formats
<b>D-Gluconic Acid</b>	K-GATE	Food additive	Rapid reaction, stable reagents
<b>D-Glucose</b>	K-GLUC K-GLUHKR/L	Common food component, very important in certain situations, e.g. diabetic products	Choice of simple formats available, based either on glucose oxidase / peroxidase, or hexokinase / G-6-PDH
<b>L-Glutamic Acid</b>	K-GLUT	Common natural food component, e.g. in cheese and tomatoes, or added as a flavouring agent, e.g. as monosodium glutamate (MSG)	Diaphorase supplied as a stabilised suspension rather than a lyophilised powder; thus less wasted enzyme
<b>Glycerol</b>	K-GCROL K-GCROLGK	Common food component, or added as a sweetener or to improve "mouth feel"	Novel tablet format offers superior stability, rapid reactions
<b>D-Lactic Acid</b>	K-DATE K-DLATE	Quality indicator of fruit and vegetable products	Rapid reaction, stable reagents
<b>L-Lactic Acid</b>	K-LATE	Quality indicator of fruit, vegetable and egg products	Rapid reaction, stable reagents. Ideal for manual and auto-analyser applications
<b>Lactose</b>	K-LACGAR K-LACSU	Common processed food component, exact amount important in "lactose free" products	Very rapid reaction for K-LACGAR (~ 5 min even at room temperature), stable reagents
<b>Maltose</b>	K-MASUG	Common food component	Rapid reaction, stable reagents
<b>Resistant Starch</b>	K-RSTAR	Starch that is not digested in the small intestine of monogastric animals	Only kit available. Rapid and robust. AOAC Method 2002.02; AACC Method 32-40.01
<b>Sucrose</b>	K-SUFRG K-SUCGL	Common food component	Choice of simple formats available, based either on glucose oxidase / peroxidase, or hexokinase / G-6-PDH
<b>Sweeteners</b>	K-ASPTM K-MANOL K-SORB	Aspartame, D-mannitol, D-sorbitol and xylitol are common sweeteners found in a variety of foods	1. K-ASPTM - novel method, only test kit available 2. K-MANOL - new method, only test kit available 3. K-SORB - diaphorase supplied as a stabilised suspension rather than a lyophilised powder; thus less wasted enzyme
<b>Total Dietary Fiber</b>	K-TDFR K-INTDF	Carbohydrate not digested in small intestine	1. K-TDFR: AOAC Methods 985.29, 991.42, 991.43 & 993.19; AACC Methods 32-05.01, 32-06.01, 32-07.01, 32-21.01, 2. K-INTDF is consistent with the CODEX Alimentarius definition of dietary fiber. AOAC Method 2009.01, 2011.25; AACC Methods 32-45.01 & 32-50.11
<b>Total Starch</b>	K-TSTA K-TSTAHK	Major food component	Rapid assay formats with options of measuring D-glucose with GOPOD reagent or with hexokinase / G-6-PDH. Stable reagents. AOAC Method 996.11; AACC Method 76-13.01; ICC Method No. 168; RACI Standard Method





Analyte	Cat. No.	Analyte Significance	Advantages of Megazyme Test Kits
<b>Acetic Acid</b>	K-ACETRM	Commonly found in feed or fermented feed	K-ACETRM is a rapid, manual assay kit employing AK and phosphotransacetylase. Stable reagents
<b>Ammonia</b>	K-AMIAR	Commonly found in feed or fermented feed	Rapid reaction rate (~ 3 min at room temperature). Ideal for manual and auto-analyser applications. Stable reagents
<b><math>\alpha</math>-Amylase</b>	K-CERA	A key enzyme in most feeds and plant products	Novel assay employing a defined oligosaccharide substrate. High sensitivity and specificity. AOAC Method 2002.01; AACC Method 22-02.01; ICC Standard Method no. 303; RACI Standard Method; CCFRA Flour Testing Working Group Method 0018
<b>Available Carbohydrates / Dietary Fiber</b>	K-ACHDF	Rapidly and slowly available sugars for digestion or fermentation	Novel procedure, stable reagents
<b>Fructan</b>	K-FRUC	Fructo-oligosaccharides in grasses and grains	Only kit available. Stable reagents; AOAC Method 999.03; AACC Method 32-32.01
<b>D-Fructose / D-Glucose</b>	K-FRUGL K-FRGLMQ K-FRGLQR	Major digestible carbohydrates in feeds	Rapid reaction times, choice of simple formats available, ideal for manual and auto-analyser applications, stable reagents
<b>Galactomannan</b>	K-GALM	Reserve carbohydrate in many legume seeds	Only kit available, stable reagents
<b><math>\beta</math>-Glucan (Barley and oats)</b>	K-BGLU	Major cell-wall polysaccharide of barley and oats	Rapid reaction, stable reagents, only enzymatic kit available. AOAC Method 995.16; AACC Method 32-23.01; EBC Methods 3.11.1, 4.16.1 and 8.11.1; ICC Standard No. 166; RACI Standard Method
<b><math>\beta</math>-Glucanase</b>	K-CELLG3	$\beta$ -Glucanase in feed	Novel assay employing a defined oligosaccharide substrate. High sensitivity, specificity and stability. Rapid reaction, ideal for manual and auto-analyser applications
<b><math>\beta</math>-Glucanase</b>	K-MBGL	Cellulase and $\beta$ -glucanase levels in feeds	Only kit available, stable reagents. RACI Standard Method
<b>L-Lactic Acid</b>	K-LATE	Commonly found in fermented feed	Rapid reaction, stable reagents. Ideal for manual and auto-analyser applications
<b>Phytic Acid</b>	K-PHYT	Found in most plant materials. Major form of bound phosphate in plant materials	Novel procedure. Rapid reaction, stable reagents
<b>Raffinose / D-Galactose</b>	K-RAFGA	Found in high levels in legume seeds. Causes discomfort and flatulence in pigs	Rapid reaction, stable reagents
<b>Resistant Starch</b>	K-RSTAR	Starch that is not digested in the small intestine of monogastric animals	Only kit available, stable reagents AOAC Method 2002.02; AACC Method 32-40.01,
<b>Total Dietary Fiber</b>	K-TDFR K-INTDF	Carbohydrate not digested in small intestine	1. K-TDFR: AOAC Methods 985.29, 991.42, 991.43 & 993.19; AACC Methods 32-05.01, 32-06.01, 32-07.01, 32-21.01, 2. K-INTDF is consistent with the CODEX Alimentarius definition of dietary fiber. AOAC Method 2009.01, 2011.25; AACC Methods 32-45.01 & 32-50.11
<b>Total Starch</b>	K-TSTA K-TSTAHK	Starch content of grain and feed	Rapid assay formats with options of measuring D-glucose with GOPOD reagent or with hexokinase / G-6-PDH. Stable reagents. AOAC Method 996.11; AACC Method 76-13.01; ICC Method No. 168; RACI Standard Method
<b>endo-<math>\beta</math>-Xylanase</b>	K-XYLS	$\beta$ -Xylanase in feed	High sensitivity, stable reagents
<b>endo-<math>\beta</math>-Xylanase</b>	S-AXBP	$\beta$ -Xylanase in feed	Sensitive, easy to use, stable reagent
<b>Protease</b>	S-AZCAS	endo-Protease added to feed	Easy to use, stable reagent



## Fermentation Industry Test Kits

Analyte	Cat. No.	Analyte Significance	Advantages of Megazyme Test Kits
<b>Acetic Acid</b>	K-ACETRM	A common fermentation product	K-ACETRM is a rapid, manual assay kit employing AK and phosphotransacetylase. Stable reagents
<b>Ammonia</b>	K-AMIAR	Commonly measured in fermentation broths	K-AMIAR has a very rapid reaction rate (~ 3 min at room temperature). Ideal for manual and auto-analyser applications. Stable reagents
<b><math>\alpha</math>-Amylase</b>	K-CERA	A major fermentation product	Novel assay employing a defined oligosaccharide substrate. High sensitivity and specificity. AOAC Method 2002.01; AACC Method 22-02.01; ICC Standard Method no. 303; RACI Standard Method; CCFA Flour Testing Working Group Method 0018
<b>L-Asparagine / L-Glutamine / Ammonia</b>	K-ASNAM	Common components of animal cell culture media	Novel product, enabling all three analytes to be determined in less than 20 min. Manual and microplate format procedures given
<b>Citric Acid</b>	K-CITR	A product of fermentation	Ideal for both manual and auto-analyser applications. Reconstituted citrate lyase stable for > 6 months at -20°C. Stable reagents
<b>Ethanol</b>	K-ETOH	Produced during alcoholic fermentation	Rapid reaction, stable reagents (AIDH supplied as a stable suspension)
<b><math>\beta</math>-Glucanase</b>	K-CELLG3	A major fermentation product	Novel assay employing a defined oligosaccharide substrate. High sensitivity, specificity and stability. Rapid reaction, ideal for manual and auto-analyser applications
<b><math>\beta</math>-Glucanase</b>	K-MBGL	A major fermentation product	Rapid reaction, stable reagents; RACI Standard Method
<b>D-Glucose</b>	K-GLUC K-GLUHK	Common component of fermentation broths	Rapid reaction, stable reagents
<b>Glucose Oxidase</b>	K-GLOX	A major fermentation product	Rapid reaction, simple format, stable reagents
<b>L-Glutamine / Ammonia</b>	K-GLNAM	Common components of animal cell culture media	Simple and rapid test kit gives values for ammonia and L-glutamine
<b>Glycerol</b>	K-GCROL K-GCROLGK	A product of fermentation	Rapid reactions, stable reagents
<b>L-Lactic Acid</b>	K-LATE	Produced predominantly from L-malic acid during malolactic fermentation	Rapid reaction, stable reagents. Ideal for manual and auto-analyser applications
<b>L-Malic Acid</b>	K-LMALR K-LMALAF K-LMALMQ K-LMALQR	Common component of fruits	All kits contain PVP to prevent tannin inhibition. 1. K-LMALR/L (manual) rapid reaction 2. K-LMALAF (auto) rapid reaction, excellent linearity 3. K-LMALMQ (manual, colorimeter based) 4. K-LMALQR (auto) liquid ready reagents
<b>Succinic Acid</b>	K-SUCC	Wine acid produced during fermentation	Rapid reaction (~ 6 min at room temperature), stable reagents
<b>Sucrose</b>	K-SUFRG K-SUCGL	Added to increase the amount of alcohol. Use only permitted in certain situations	Choice of simple formats available, based either on glucose oxidase / peroxidase, or hexokinase / G-6-PDH
<b>Urea</b>	K-URAMR	Source of Yeast Available Nitrogen (YAN) and precursor of the carcinogen ethyl carbamate. Over-supplementation with diammonium phosphate (DAP) can result in elevated levels	Simple, very rapid (both urea and ammonia measured in < 10 min at room temperature) and sequential / efficient (only one cuvette required per sample)
<b><math>\alpha</math>-Amylase</b>	T-AMZ200	A product of fermentation	Rapid reaction, stable reagent AACC Method 22.05; RACI Standard Method
<b>endo-Arabinanase</b>	T-ARZ200	A product of fermentation	Rapid reaction, stable reagent
<b><math>\beta</math>-Glucanase</b>	S-ABG100	A product of fermentation	Rapid reaction, stable reagent
<b>Pullulanase</b>	S-RPUL	A product of fermentation	Rapid reaction, stable reagent
<b>endo-<math>\beta</math>-Xylanase</b>	S-AXBP	A product of fermentation	Rapid reaction, stable reagent



## Wine Industry Test Kits



Analyte	Cat. No.	Analyte Significance	Advantages of Megazyme Test Kits
<b>Acetaldehyde</b>	K-ACHYD	A sensory compound that adds flavour and complexity, but spoils wine at high concentrations	AcDH supplied as a stabilised solution rather than a lyophilised powder; thus less wasted enzyme
<b>Acetic Acid</b>	K-ACET K-ACETAF K-ACETAK K-ACETRM K-ACETGK	A sensory compound that adds flavour and complexity in small amounts, but spoils wine at high concentrations. Produced naturally by yeast in small amounts and by spoilage organisms such as <i>Acetobacter aceti</i> in large quantities. This is the predominant of the acids comprising ~ 85% volatile acidity (VA)	All kits contain PVP to prevent tannin inhibition. K-ACET (manual, efficient) contains stable ACS suspension. K-ACETAF (auto) used to prepare very stable R1 and R2. K-ACETAK (auto) / K-ACETRM (manual) are very rapid acetate kinase (AK) based kits with excellent linearity. K-ACETGK is a new rapid, auto-analyser assay kit employing AK and phosphotransacetylase. Stable reagents
<b>Ammonia</b>	K-AMIAR K-LARGE	Most important inorganic source of Yeast Available Nitrogen (YAN)	Novel enzyme employed is not inhibited by tannins, endpoint reaction time ~ 3 min. Ideal for manual and auto-analyser applications
<b>L-Arginine</b>	K-LARGE	Most important amino acid in grape juice with respect to YAN	Simple and rapid test kit gives sequential values for ammonia, urea and L-arginine. No tannin inhibition
<b>L-Ascorbic Acid</b>	K-ASCO	Present naturally in grapes and can be added as an anti-oxidant	Rapid reaction, stable reagents
<b>Citric Acid</b>	K-CITR	Naturally present in small amounts; large amounts indicate addition for acidification (EU limit is 1 g/L)	Ideal for both manual and auto-analyser applications. Reconstituted citrate lyase stable for > 6 months at -20°C. Stable reagents
<b>Ethanol</b>	K-ETOH	Produced during alcoholic fermentation. Amounts > 17.5% (v/v) indicate supplementation	Rapid reaction, stable reagents (AIDH supplied as a stable suspension)
<b>D-Fructose / D-Glucose</b>	K-FRUGL K-FRGLMQ K-FRGLQR	Grape quality indicator. One of the two principle fermentable sugars of grape juice	Contains PVP to prevent tannin inhibition. Ideal for manual and auto-analyser use. Stable reagents
<b>D-Gluconic Acid</b>	K-GATE	Grape quality indicator for the production of certain wines	Rapid reaction, stable reagents
<b>Glycerol</b>	K-GCROL K-GCROLGK	Quality indicator of finished wine, important for "mouth feel"	Novel tablet format offers superior stability, rapid reaction
<b>D-Lactic Acid</b>	K-DATE K-DLATE	Produced predominantly by lactic acid spoilage bacteria	Rapid reaction, stable reagents
<b>L-Lactic Acid</b>	K-LATE K-DLATE	Produced predominantly from L-malic acid during malolactic fermentation	Rapid reaction, stable reagents. Ideal for manual and auto-analyser applications
<b>D-Malic Acid</b>	K-DMAL	Only present in significant quantities in adulterated wine	D-MDH supplied as a stabilised suspension rather than a lyophilised powder; thus less wasted enzyme
<b>L-Malic Acid</b>	K-LMALR K-LMALAF K-LMALMQ K-LMALQR	Grape quality indicator. Very important grape acid, converted to less acidic L-lactic acid during malolactic fermentation	All kits contain PVP to prevent tannin inhibition. 1. K-LMALR/L (manual) rapid reaction 2. K-LMALAF (auto) rapid reaction, excellent linearity 3. K-LMALMQ (manual, colorimeter based) 4. K-LMALQR (auto) liquid ready reagent
<b>Primary Amino Nitrogen (NOPA)</b>	K-PANOPA	Primary amino nitrogen (PAN) is the most important organic source of YAN	Novel kit, rapid reaction, stable reagents, simple format
<b>D-Sorbitol</b>	K-SORB	High levels indicate addition of fruit	Diaphorase supplied as a stabilised suspension rather than a lyophilised powder; thus less wasted enzyme
<b>Succinic Acid</b>	K-SUCC	Wine acid produced during fermentation	Rapid reaction (~ 6 min even at room temperature), stable reagents
<b>Sucrose</b>	K-SUFRG K-SUCGL	Added to increase the amount of alcohol. Use only permitted in certain situations	Choice of simple formats available, based either on glucose oxidase / peroxidase, or hexokinase / G-6-PDH
<b>Sulphite</b>	K-SULPH K-TSULPH K-ETSULPH	Sulphites are used as an essential additive in the control of microbial contamination during aging and also to protect the wine against detrimental "oxidative and enzymatic browning"	Choice of simple formats available, based either on liquid ready reagent chemical reactions (K-SULPH & K-TSULPH) or an enzymatic reaction (K-ETSULPH). Stable reagents
<b>Tartaric Acid</b>	K-TART	Occurs naturally in grapes and is one of the most prevalent organic acids. Key indicator of total (titratable) acidity (TA)	Stable liquid ready reagents. Simple, rapid chemical reaction for manual, auto-analyser and microplate formats
<b>Urea</b>	K-URAMR	Source of YAN and precursor of the carcinogen ethyl carbamate. Over-supplementation with DAP can result in elevated levels	Simple, very rapid (both urea and ammonia measured in < 10 min at room temperature) and sequential / efficient (only one cuvette required per sample)





## Brewing Industry Test Kits



Analyte	Cat. No.	Analyte Significance	Advantages of Megazyme Test Kits
$\alpha$ -Amylase	K-CERA	A key indicator of malt quality	Novel assay employing a defined oligosaccharide substrate. High sensitivity and specificity. AOAC Method 2002.01; AACC Method 22-02.01; ICC Standard Method no. 303; RACI Standard Method; CCFRA Flour Testing Working Group Method 0018
$\beta$ -Amylase	K-BETA3	A key indicator of malt quality	Only kit available. Stable reagents; RACI Standard Method
$\beta$ -Glucan (Barley and oats)	K-BGLU	Major cell-wall polysaccharide of barley and oats	Rapid reaction, stable reagents, only enzymatic kit available. AOAC Method 995.16; AACC Method 32-23.01; EBC Methods 3.11.1, 4.16.1 and 8.11.1; ICC Standard Method No. 166; RACI Standard Method
$\beta$ -Glucanase	K-MBGL	$\beta$ -Glucanase level in malt	Rapid reaction, stable reagents, only enzymatic kit available; RACI Standard Method
D-Glucose	K-GLUC K-GLUHKR/L	Major component of fermentation mixture	Rapid reaction, stable reagents
Malt Amylase	K-MALTA	Measurement of $\alpha$ -/ $\beta$ -amylase. Key indicators of malt quality	Combination of both K-CERA and K-BETA3
Total Starch	K-TSTA K-TSTAHK	Starch content of grain and feed	Rapid assay formats with options of measuring D-glucose with GOPOD reagent or with hexokinase / G-6-PDH. Stable reagents. AOAC Method 996.11; AACC Method 76-13.01; ICC Method No. 168; RACI Standard Method
Alpha-Amylase	T-AMZBG200	Allows measurement of $\alpha$ -amylase in pre-harvest sprouted barley	Novel procedure. Rapid reaction, stable reagent
$\beta$ -Glucanase	T-BGZ200	Key enzyme in hydrolysis of malt $\beta$ -glucans	Novel substrate. Rapid reaction, stable reagent; RACI Standard Method
Limit-Dextrinase	T-LDZ200	Key enzyme in hydrolysis of 1,6-linkages in starch and branched malto-dextrins	Novel substrate. Rapid reaction, stable reagent; RACI Standard Method
endo- $\beta$ -Xylanase	T-XAX200	Key enzyme in hydrolysis of malt xylans	Novel substrate. Rapid reaction, stable reagent





Analyte	Cat. No.	Analyte Significance	Advantages of Megazyme Test Kits
<b>Acetaldehyde</b>	K-ACHYD	One of the most important aroma compounds in yogurt, responsible for the characteristic taste. Also formed in milk during storage	AIDH supplied as a stabilised solution rather than a lyophilised powder, thus less wasted enzyme. Stable reagents
<b>Acetic Acid</b>	K-ACET K-ACETAF K-ACETAK K-ACETRM K-ACETGK	Fermentation product of yogurt and cheese	All kits contain PVP to prevent tannin inhibition. K-ACET (manual, efficient) contains stable ACS suspension. K-ACETAF (auto) used to prepare very stable R1 and R2. K-ACETAK (auto) / K-ACETRM (manual) are very rapid acetate kinase (AK) based kits with excellent linearity. K-ACETGK is a new rapid, auto-analyser assay kit employing AK and phosphotransacetylase. Stable reagents
<b>Ammonia</b>	K-AMIA K-AMIAR	Important indicator of the hygienic quality (microbial load) of milk	K-AMIAR has a very rapid reaction rate (~ 3 min at room temperature). Ideal for manual and auto-analyser applications, stable reagents
<b>L-Ascorbic Acid</b>	K-ASCO	Antioxidant present in dairy products. Permitted additive	Rapid reaction, stable reagents
<b>Aspartame</b>	K-ASPTM	Common milkshake and yogurt sweetener	Rapid reaction, stable reagents, only enzymatic kit available
<b>Citric Acid</b>	K-CITR	Important quality indicator of milk, especially for butter and cheese production. Permitted additive	Ideal for both manual and auto-analyser applications. Reconstituted citrate lyase stable for > 6 months at -20°C, stable reagents
<b>Ethanol</b>	K-ETOH	Produced during the fermentation of kefir	Rapid reaction, stable reagents (AIDH supplied as a stable suspension)
<b>Formic Acid</b>	K-FORM	Minor acid in dairy products	FDH supplied as a stabilised suspension rather than a lyophilised powder, thus less wasted enzyme, stable reagents
<b>D-Fructose / D-Glucose</b>	K-FRUGL K-FRGLMQ	Common milkshake and yogurt sweetener	Rapid reaction times, choice of simple formats available, ideal for manual and auto-analyser applications, stable reagents
<b>D-Gluconic Acid</b>	K-GATE	Weak organic acid found in dairy products. High levels found in certain cheeses	Rapid reaction, stable reagents
<b>D-Glucose</b>	K-GLUC K-GLUHKR/L	Low levels expected in unprocessed / unadulterated milk and in cheese. Useful marker when producing lactose depleted dairy products	Choice of simple formats available, based either on glucose oxidase / peroxidase, or hexokinase / G-6-PDH, stable reagents
<b>L-Glutamic Acid</b>	K-GLUT	Found in high concentrations, especially in cheese	No wasted diaphorase solution (stable suspension supplied), stable reagents
<b>D-Lactic Acid</b>	K-DATE	Quality indicator of milk, yogurt and cheese	Rapid reaction, stable reagents
<b>L-Lactic Acid</b>	K-LATE	Quality indicator of fresh milk. High levels in yogurt and cheese	Rapid reaction, stable reagents. Ideal for manual and auto-analyser applications
<b>D-/L-Lactic Acid</b>	K-DLATE	Quality indicator of fresh milk, yogurt and cheese	Rapid reaction, flexible concurrent format, stable reagents
<b>Lactose / D-Galactose</b>	K-LACGAR	Key quality (value) indicator of milk	Very rapid reaction (~ 5 min even at room temperature), stable reagents
<b>D-Sorbitol / Xylitol</b>	K-SORB	Dairy product sweetener	No wasted diaphorase solution (stable suspension supplied), stable reagents
<b>Succinic Acid</b>	K-SUCC	Minor dairy acid	Rapid reaction (~ 6 min even at room temperature), stable reagents
<b>Sucrose</b>	K-SUFRG K-SUCGL	Not present naturally in dairy products	Choice of simple formats available, based either on glucose oxidase / peroxidase, or hexokinase / G-6-PDH, stable reagents
<b>Urea</b>	K-URAMR	Quality indicator of milk, especially that used for cheese production. Used as a metabolic marker of bovine blood urea levels	Simple, very rapid (both urea and ammonia measured in < 10 min at room temperature) and sequential / efficient (only one cuvette required per sample)





# OVERVIEW OF ASSAY METHODS







## Acetaldehyde

Cat. No. **K-ACHYD**

**UV-method for the determination of Acetaldehyde in foodstuffs, beverages and other materials**

**Principle:**

(aldehyde dehydrogenase)



<b>Kit size:</b>	50 assays (manual) / 500 (microplate) / 500 (auto-analyser)
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Reaction time:</b>	~ 4 min
<b>Detection limit:</b>	0.18 mg/L
<b>Application examples:</b>	Wine, champagne, beer, liqueurs, brandy, dairy products (e.g. yogurt), bread, fruit juices, soft drinks, cocoa, vegetable and fruit products, coffee, and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	Methods based on this principle have been accepted by MEBAK

**Advantages**

- No wasted aldehyde dehydrogenase solution (stable suspension supplied)
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Simple format
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats



## Acetic Acid (Acetyl-CoA synthetase analyser format)

Cat. No. **K-ACETAF**

**Analyser format UV-method for the determination of Acetic Acid in foodstuffs, beverages and other materials**

**Principle:**


(acetyl-CoA synthetase)



(citrate synthase)



(L-malate dehydrogenase)



<b>Kit size:</b>	141.6 mL of prepared reagent (R1 + R2)
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Reaction time:</b>	~ 15 min
<b>Detection limit:</b>	10 mg/L (recommended assay format)
<b>Application examples:</b>	Wine, beer, fruit and fruit juices, soft drinks, vinegar, vegetables, pickles, dairy products (e.g. cheese), meat, fish, bread, bakery products (and baking agents), ketchup, soy sauce, mayonnaise, dressings, paper (and cardboard), tea, pharmaceuticals (e.g. infusion solutions), feed and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	Methods based on this principle have been accepted by EN, ISO, ICUMSA, IFU, MEBAK

**Advantages**

- No wasted ACS solution (stable suspension supplied)
- PVP incorporated to prevent tannin inhibition
- Very stable reagent when prepared for auto-analyser applications (> 3 days at 4°C)
- Linear calibration up to 30 µg/mL of acetic acid in final reaction solution
- Validated by the University of Wine, Suze la Rousse, France
- Very competitive price (cost per mL of reagent)
- All reagents stable for > 2 years after preparation



## Acetic Acid (Acetyl-CoA synthetase manual format)

Cat. No. **K-ACET**

**Manual format UV-method for the determination of Acetic Acid in foodstuffs, beverages and other materials**

**Principle:**


(acetyl-CoA synthetase)



(citrate synthase)



(L-malate dehydrogenase)



**Kit size:** 53 assays

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 14 min

**Detection limit:** 0.14 mg/L

**Application examples:** Wine, beer, fruit and fruit juices, soft drinks, vinegar, vegetables, pickles, dairy products (e.g. cheese), meat, fish, bread, bakery products (and baking agents), ketchup, soy sauce, mayonnaise, dressings, paper (and cardboard), tea, pharmaceuticals (e.g. infusion solutions), feed, and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by EN, ISO, ICUMSA, IFU, MEBAK

### Advantages

- No wasted ACS solution (stable suspension supplied)
- PVP incorporated to prevent tannin inhibition
- All reagents stable for > 2 years after preparation
- Very competitive price (cost per test)
- *Mega-Cal*™ software tool is available from our website for hassle-free raw data processing




## Acetic Acid (Acetate kinase analyser format)

Cat. No. **K-ACETAK**


**Analyser format UV-method for the determination of Acetic Acid in foodstuffs, beverages and other materials**

**Principle:**

(acetate kinase)



(pyruvate kinase)



(D-lactate dehydrogenase)



**Kit size:** 550 assays

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 10 min

**Detection limit:** 10 mg/L (recommended assay format)

**Application examples:** Wine, beer, fruit and fruit juices, soft drinks, vinegar, vegetables, pickles, dairy products (e.g. cheese), meat, fish, bread, bakery products (and baking agents), ketchup, soy sauce, mayonnaise, dressings, paper (and cardboard), tea, pharmaceuticals (e.g. infusion solutions), feed, and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** *Improved method*

### Advantages

- Very stable reagent when prepared for auto-analyser applications (> 7 days at 4°C)
- PVP incorporated to prevent tannin inhibition
- Linear calibration ( $R^2 \sim 0.9995$ ) up to 30 µg/mL of acetic acid in final reaction solution
- Validated by the University of Wine, Suze la Rousse, France
- Very rapid reaction
- Very competitive price (cost per mL of reagent)
- All reagents stable for > 2 years
- Extended cofactors stability




## Acetic Acid (Analyser format)

Cat. No. **K-ACETGK**

**Analyser format UV-method for the determination of Acetic Acid in foodstuffs, beverages and other materials**

### Principle:

(acetate kinase)




(PTA)



(ADP-GK)



(G6P-DH)



**Kit size:**

500 assays

**Method:**

Spectrophotometric at 340 nm

**Reaction time:**

8 min at 25°C or 5 min at 37°C

**Detection limit:**

1.8 g/L (recommended assay format)

**Application examples:**

Wine, beer, fruit and fruit juices, soft drinks, vinegar, vegetables, pickles, dairy products (e.g. cheese), meat, fish, bread, bakery products (and baking agents), ketchup, soy sauce, mayonnaise, dressings, paper (and cardboard), tea, pharmaceuticals (e.g. infusion solutions), feed, and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:**

**Improved method**




## Acetic Acid (Rapid manual format)

Cat. No. **K-ACETRM**

**Manual format UV-method for the determination of Acetic Acid in foodstuffs, beverages and other materials**

### Principle:


(acetate kinase)




(phosphotransacetylase)



(pyruvate kinase)



(D-lactate dehydrogenase)



**Kit size:**

72 assays (manual) / 720 (microplate)

**Method:**

Spectrophotometric at 340 nm

**Reaction time:**

~ 4 min

**Detection limit:**

0.063 mg/L

**Application examples:**

Wine, beer, fruit and fruit juices, soft drinks, vinegar, vegetables, pickles, dairy products (e.g. cheese), meat, fish, bread, bakery products (and baking agents), ketchup, soy sauce, mayonnaise, dressings, paper (and cardboard), tea, pharmaceuticals (e.g. infusion solutions), feed and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:**

**Improved method**

### Advantages

- Excellent reagent stability
- > 7 days at 4°C or > 2 years at -20°C when prepared for auto-analyser applications
- > 2 years as supplied
- Very rapid reaction (~ 5 min at 37°C)
- Linear calibration (R2 ~ 0.997 up to 1.8 g/L sample)

### Advantages

- Improved assay format (only two absorbance readings required)
- All reagents stable for > 2 years after preparation
- PVP incorporated to prevent tannin inhibition
- Very rapid reaction (~ 4 min)
- Mega-Calc™ software tool is available from our website for hassle-free raw data processing
- Very competitive price (cost per test)
- Suitable for Manual and Microplate formats



Food



Feed



Fermentation



Wine



Dairy


## Ammonia (Rapid)

Cat. No. K-AMIAR

**UV-method for the determination of Ammonia in foodstuffs, beverages and other materials**

### Principle:

(microbial glutamate dehydrogenase)



**Kit size:** 96 assays (manual) / 960 (microplate) / 960 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 3 min

**Detection limit:** 0.07 mg/L

**Application examples:** Grape juice, wine, fruit juices, soft drinks, dairy products (e.g. milk), dietetic food, soy sauce, eggs and egg products, cheese, meat, processed meat, seafood, bakery products (and baking agents), fertilisers, pharmaceuticals, tobacco, cosmetics, water, Kjeldahl analysis, paper (and cardboard), water and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by MEBAK

### Advantages

- Very rapid reaction due to use of uninhibited glutamate dehydrogenase
- Enzyme supplied as stabilised suspension
- Very competitive price (cost per test)
- All reagents stable for > 2 years as supplied
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats



Food



Feed



Fermentation



Brewing



Biofuels


## α-Amylase ("Ceralpha" Method)

Cat. No. K-CERA

**Colourimetric method for the determination of Alpha-Amylase in foodstuffs, feed and fermentation products**

### Principle:

(α-amylase)



(thermostable α-glucosidase)



(alkaline solution)



Note: PNP = 4-nitrophenol

**Kit size:** 100 / 200 assays

**Method:** Spectrophotometric at 400 nm

**Total assay time:** ~ 20 min

**Detection limit:** 0.05 U/mL

**Application examples:** Cereal flours, fermentation broths and other materials

**Method recognition:** AOAC (Method 2002.01), AACC (Method 22.02.01), ICC (Standard No. 303), RACI (Standard Method), and CCFRA (Flour Testing Working Group Method 0018)

### Advantages

- Very cost effective
- All reagents stable for > 2 years after preparation
- Very specific
- Simple format
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standard included





Food



Feed



Fermentation



Brewing



Biofuels

## $\alpha$ -Amylase ("Sprout Damaged" Method)

Cat. No. **K-AMYLS**

Highly sensitive colourimetric method for the determination of  $\alpha$ -Amylase in sprout damaged grain

### Principle:

- ( $\alpha$ -amylase)
- (1) Ethylidene- $G_7$ - $\alpha$ -PNP +  $H_2O \rightarrow$  Ethylidene- $G_x$  +  $G_{(7-x)}$ - $\alpha$ -PNP
- (thermostable  $\alpha$ -glucosidase)
- (2)  $G_{(7-x)}$ - $\alpha$ -PNP +  $H_2O \rightarrow$  D-glucose + PNP
- (alkaline solution)
- (3) PNP  $\rightarrow$  phenolate ion (yellow colour)

Note: PNP = 4-nitrophenol

<b>Kit size:</b>	80 / 160 assays (manual) / 640 (auto-analyser)
<b>Method:</b>	Spectrophotometric at 400 nm
<b>Total assay time:</b>	~ 5 min
<b>Detection limit:</b>	0.05 U/mL
<b>Application examples:</b>	Sprout damaged wheat grain
<b>Method recognition:</b>	Novel method

### Advantages

- Extremely high sensitivity  
- 2.4-fold increase over Ceralpha (K-CERA)
- Very cost effective
- All reagents stable for > 2 years after preparation
- Very specific
- Simple format
- Mega-Calc™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual and auto-analyser formats



Food



Feed



Fermentation



Brewing

## $\beta$ -Amylase ("Betamyl-3" Method)

Cat. No. **K-BETA3**

Colourimetric method for the determination of  $\beta$ -Amylase in cereal grains, malt, food, beverages and fermentation products

### Principle:

- ( $\beta$ -amylase)
- (1)  $G_3$ - $\beta$ -PNP +  $H_2O \rightarrow$   $G_2$  +  $G$ - $\beta$ -PNP
- ( $\beta$ -glucosidase)
- (2)  $G$ - $\beta$ -PNP +  $H_2O \rightarrow$  D-glucose + PNP
- (alkaline solution)
- (3) p-Nitrophenol  $\rightarrow$  phenolate ion (yellow colour)

Note: PNP = 4-nitrophenol

<b>Kit size:</b>	100 / 200 assays
<b>Method:</b>	Spectrophotometric at 400 nm
<b>Reaction time:</b>	~ 10 min
<b>Detection limit:</b>	0.05 U/mL of sample solution
<b>Application examples:</b>	Cereal flours, malts and other materials
<b>Method recognition:</b>	Modification of RACI (Standard Method)

### Advantages

- Very cost effective
- All reagents stable for > 2 years as supplied
- Only enzymatic kit available
- Very specific
- Simple format
- Rapid reaction
- Mega-Calc™ software tool is available from our website for hassle-free raw data processing
- Standard included



Food







Feed

## Amylose / Amylopectin

Cat. No. **K-AMYL**

Colourimetric method for the determination of Amylose and Amylopectin in cereals, food and feed

### Principle:



<b>Kit size:</b>	100 assays
<b>Method:</b>	Spectrophotometric at 510 nm
<b>Total assay time:</b>	~ 120 min
<b>Detection limit:</b>	Amylose 5-95% of total starch content
<b>Application examples:</b>	Cereal starches, flours, pure starches and foods
<b>Method recognition:</b>	<b>Novel method</b>

### Advantages

- Very cost effective (cost per test)
- All reagents stable for > 12 months after preparation
- Only enzymatic kit available
- Accurate and reliable amylose / amylopectin ratio determination
- Simple format
- Standard included



Food






Feed

## Amyloglucosidase / Glucoamylase

Cat. No. **K-AMG**

Colourimetric method for the determination of Amyloglucosidase

### Principle:



Note: PNP = 4-nitrophenol

<b>Kit size:</b>	150 / 300 assays
<b>Method:</b>	Spectrophotometric at 400 nm
<b>Total assay time:</b>	~ 10 min
<b>Detection limit:</b>	0.05 U/mL
<b>Application examples:</b>	Industrial enzyme preparations
<b>Method recognition:</b>	<b>Novel method</b>

### Advantages

- Highly sensitive
- Very cost effective
- All reagents stable for > 2 years after preparation
- Simple format
- *Mega-Cal<sup>TM</sup>* software tool is available from our website for hassle-free raw data processing
- Standard included



Feed



Biofuels


## Arabinan

Cat. No. K-ARAB

UV-method for the determination of Arabinan in plant materials and juices

### Principle:

(endo-arabinanase +  $\alpha$ -L-arabinofuranosidase)



(galactose mutarotase)



( $\beta$ -galactose dehydrogenase)



<b>Kit size:</b>	100 assays
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Reaction time:</b>	~ 10 min
<b>Detection limit:</b>	1.3 mg/L
<b>Application examples:</b>	Fruit juices and other materials
<b>Method recognition:</b>	<b>Novel method</b>

### Advantages

- Very rapid reaction due to inclusion of galactose mutarotase (patented technology)
- Very cost effective
- All reagents stable for > 2 years after preparation
- Only enzymatic kit available
- Very specific
- Simple format
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included



Fermentation



Biofuels

## L-Arabinose / D-Galactose (Rapid)

Cat. No. K-ARGA

UV-method for the determination of L-Arabinose and D-Galactose in hydrolysed plant products

### Principle:

(galactose mutarotase)



( $\beta$ -galactose dehydrogenase)



( $\beta$ -galactose dehydrogenase)



<b>Kit size:</b>	115 assays (manual) / 1150 (microplate) / 1150 (auto-analyser)
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Reaction time:</b>	~ 5 min
<b>Detection limit:</b>	1.3 mg/L
<b>Application examples:</b>	Analysis of hydrolysates of oligo- and polysaccharides (e.g. arabinan, arabinoxylan, galactan, arabinogalactan), milk, dairy products, foods containing milk (e.g. dietetic foods, bakery products, baby food, chocolate, sweets and ice-cream), food additives (e.g. sweeteners), cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	<b>Novel method</b>

### Advantages

- Very rapid reaction due to inclusion of galactose mutarotase (patented technology)
- Very cost effective
- All reagents stable for > 2 years after preparation
- Only enzymatic kit available
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



Fermentation



Wine


## L-Arginine / Urea / Ammonia (Rapid)

Cat. No. **K-LARGE**

**UV-method for the determination of L-Arginine, Urea and Ammonia in grape juice, must and wine**

**Principle:**


(microbial glutamate dehydrogenase)



(urease)



(arginase)



**Kit size:** 50 assays of each

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 20 min [ammonia (2 min), urea (6 min), L-arginine (7 min)]

**Detection limit:** 0.07 mg/L (ammonia), 0.13 mg/L (urea), 0.37 mg/L (L-arginine)

**Application examples:** Grape juice, wine must, wine and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** *Improved method*

**Advantages**

- Improved assay format
- Very rapid reactions due to use of uninhibited glutamate dehydrogenase
- All enzymes supplied as stabilised suspensions
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability



Food



Fermentation



Wine

## L-Ascorbic Acid

Cat. No. **K-ASCO**


**Colourimetric method for the determination of L-Ascorbic Acid in foodstuffs, feed, wine and other materials**

**Principle:**

(5-methylphenazinium methosulphate)



(ascorbic acid oxidase)



**Kit size:** 40 assays (manual) / 400 (microplate) / 400 (auto-analyser)

**Method:** Spectrophotometric at 578 nm

**Reaction time:** ~ 8 min

**Detection limit:** 0.175 mg/L

**Application examples:** Wine, beer, fruit juices, soft drinks, jam, milk, dairy products (e.g. cheese), dietetic foods, baby foods, processed meat, baking additives, fruit and vegetables (e.g. tomato and potato), pharmaceuticals, feed and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by MEBAK

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 6 months after preparation
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats





Food




Fermentation

## L-Asparagine / L-Glutamine / Ammonia (Rapid)

Cat. No. **K-ASNAM**

**UV-method for the determination of L-Asparagine, L-Glutamine and Ammonia in potatoes, foodstuffs and cell culture media**

**Principle:**



**Kit size:** 50 assays of each

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 20 min

**Detection limit:** 0.50 mg/L (L-asparagine)  
0.54 mg/L (L-glutamine)  
0.06 mg/L (ammonia)

**Application examples:** Potatoes, potato products, vegetables, cereals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** *Novel method*

### Advantages

- Very rapid reaction due to use of uninhibited glutamate dehydrogenase
- All enzymes supplied as stabilised suspensions
- Only kit available
- Very cost effective
- All reagents stable for > 2 years after preparation
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included







Food

## Aspartame

Cat. No. **K-ASPTM**

**UV-method for the determination of Aspartame (and breakdown products) in foodstuffs, beverages and other materials**

**Principle:**



**Kit size:** 50 assays (manual) / 500 (microplate) / 500 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 5 min

**Detection limit:** 0.57 mg/L

**Application examples:** Soft drinks, artificial sweeteners, candies, mints, chewing gum, dietetic products, jam, chocolate and other materials

**Method recognition:** *Novel method*

### Advantages

- Very cost effective
- All reagents stable for > 12 months after preparation
- Only enzymatic kit available
- Measures aspartame and breakdown products (L-aspartate and aspartame acid)
- Very specific
- Very rapid reaction
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



## Available Carbohydrates / Dietary Fiber

Cat. No. **K-ACHDF**


An integrated procedure for the measurement of **Available Carbohydrates** and **Dietary Fiber** in cereal products, fruit and vegetables and food products

### Principle (Dietary Fiber):

( $\alpha$ -amylase + amyloglucosidase)




(protease)



(3) **Dietary fiber determined gravimetrically following alcohol precipitation**

### Principle (Available Carbohydrates):

(sucrase / maltase +  $\beta$ -galactosidase)



(PGI, hexokinase and glucose-6-phosphate dehydrogenase)



**Kit size:**

100 assays of each

**Application examples:**

Food ingredients, food products and other materials

**Method recognition:**

Dietary Fibre - AOAC (Methods 985.29, 991.42, 991.43 and 993.19) and AACC (Methods 32-05.01, 32-07.01 and 32-21.01)

### Advantages

- Very cost effective
- All reagents stable for > 2 years after preparation
- High purity / standardised enzymes employed
- Only kit available
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Simple format




## endo-1,4- $\beta$ -Glucanase (cellulase)

Cat. No. **K-CELLG3**

Colourimetric method for the determination of **endo-1,4- $\beta$ -glucanase (cellulase)** in enzyme preparations and fermentation products

### Principle:

(endo-1,4- $\beta$ -glucanase)



(thermostable  $\beta$ -glucosidase)



(alkaline solution)



Note: CNP = 2-Chloro-4-nitrophenol

**Kit size:**

180 / 360 assays

**Method:**

Spectrophotometric at 400 nm

**Total assay time:**

~ 20 min

**Detection limit:**

0.05 U/mL

**Application examples:**

Fermentation broths, industrial enzyme preparations, biofuels research

**Method recognition:**

**Novel method**

**Note:**

The analogous fluorimetric reagent (R-CELLF) is also available with 10-fold greater sensitivity.

### Advantages

- Cost effective
- All reagents stable for > 2 years after preparation
- Completely specific for cellulase (endo-1,4- $\beta$ -glucanase). The substrate is not hydrolysed by  $\beta$ -glucosidase, cellobiohydrolase or any other enzymes tested
- Kinetic assays possible due to significant phenolate ion presence (and UV absorbance) at pH 5-6
- Simple format. Well suited to automation
- Standard included



Food



Fermentation



Wine

## Citric Acid

Cat. No. K-CITR

**UV-method for the determination of Citric Acid in foods, beverages and other materials**

**Principle:**

(citrate lyase)

**(1) Citrate → oxaloacetate + acetate**

(L-malate dehydrogenase)

**(2) Oxaloacetate + NADH + H<sup>+</sup> → L-malate + NAD<sup>+</sup>**

(D-lactate dehydrogenase)

**(3) Pyruvate + NADH + H<sup>+</sup> → D-lactate + NAD<sup>+</sup>**

**Kit size:** 72 assays (manual) / 720 (microplate) / 840 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 5 min

**Detection limit:** 0.921 mg/L

**Application examples:** Grape juice, wine, beer, fruit juices, soft drinks, tea, dairy products (e.g. cheese), meat, processed meat, vegetable and fruit products, bakery products, paper, pharmaceuticals, cosmetics and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by MEBAK, OIV, EU, ISO2963, AOAC and IFU22

(**Note:** If the enzyme oxaloacetate decarboxylase is present in the sample, some of the oxaloacetate product is converted to pyruvate. Therefore, to ensure citric acid is measured quantitatively, D-lactate dehydrogenase (D-LDH) is employed to efficiently convert any pyruvate produced into D-lactate and NAD<sup>+</sup>).

**Advantages**

- Reconstituted citrate lyase stable for 4 weeks at 4°C / 6 months at -20°C
- Buffer / cofactor / enzyme tablets for efficient use of kit components
- PVP incorporated to prevent tannin inhibition
- Very competitive price (cost per test)
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats



Food



Fermentation



Wine



Brewing



Biofuels

## Ethanol

Cat. No. K-ETOH

**UV-method for the determination of Ethanol in foodstuffs, beverages, and other materials**

**Principle:**

(alcohol dehydrogenase)

**(1) Ethanol + NAD<sup>+</sup> ⇌ acetaldehyde + NADH + H<sup>+</sup>**

(aldehyde dehydrogenase)

**(2) Acetaldehyde + NAD<sup>+</sup> + H<sub>2</sub>O → acetic acid + NADH + H<sup>+</sup>**

**Kit size:** 60 assays (manual) / 600 (microplate) / 600 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 5 min

**Detection limit:** 0.093 mg/L

**Application examples:** Wine, beer, cider, alcoholic fruit juices, spirits, liqueurs, low-alcoholic / non-alcoholic beverages, pickles, fruit and fruit juice, chocolate products, vinegar, jam, bread and bakery products, honey, soy sauce, dairy products, cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by IFU, EBC, MEBAK, ASBC

**Advantages**

- Simple format – aldehyde dehydrogenase supplied as stable suspension
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Rapid reaction
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats



## Total Dietary Fiber

Cat. No. **K-TDFR**

For the determination of Total Dietary Fiber in cereal products, foodstuffs, feeds and other materials

### Principle:

( $\alpha$ -amylase + amyloglucosidase)

(1) Starch + H<sub>2</sub>O → D-glucose

(protease)

(2) Protein + H<sub>2</sub>O → peptides

(3) Dietary fiber determined gravimetrically following alcohol precipitation

(4) Ash and residual protein determined on DF residues and subtracted

<b>Kit size:</b>	200 assays
<b>Method:</b>	Hydrolysis / removal of non-dietary fibre components
<b>Total assay time:</b>	~ 100 min
<b>Detection limit:</b>	0.5-100% of sample weight
<b>Application examples:</b>	Food ingredients, food products and other materials
<b>Method recognition:</b>	AOAC (Methods 985.29, 991.42, 991.43 and 993.19), AACC (Methods 32-05.01, 32-06.01, 32-07.01 and 32-21.01) and CODEX (Type I Method)

### Advantages

- Very competitive price (cost per test)
- All reagents stable for > 2 years
- High purity / standardised enzymes employed
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Simple format



## Total Dietary Fiber (Integrated)

Cat. No. **K-INTDF**

For the determination of Total Dietary Fiber in cereal products, foodstuffs, feeds and other materials. **CODEX Type I method (2011)**

### Principle:

(Pancreatic  $\alpha$ -amylase + amyloglucosidase)

(1) Non-resistant starch + H<sub>2</sub>O → D-glucose

(protease)

(2) Protein + H<sub>2</sub>O → peptides

(3) IDF (including resistant starch) and alcohol precipitated soluble DF (SDFP) determined gravimetrically

(4) Alcohol soluble DF (SDFS) determined by HPLC

(5) Ash and residual protein determined on DF residues and subtracted

<b>Kit size:</b>	100 assays
<b>Method:</b>	Hydrolysis / removal of non-dietary fibre components
<b>Total assay time:</b>	~ 3 hr work (over 2 days)
<b>Detection limit:</b>	0.5-100% of sample weight
<b>Application examples:</b>	Food ingredients, food products and other materials
<b>Method recognition:</b>	AOAC (Methods 2009.01; 2011.25) AACC (Method 32-45.01; 32-50.01) CODEX (Type I method)

### Advantages

- The only method that is consistent with the CODEX Alimentarius definition of dietary fibre
- High purity / standardised enzymes employed
- All reagents stable for > 2 years
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Very competitive price (cost per test)





## Formic Acid

Cat. No. **K-FORM**

**UV-method for the determination of Formic Acid in foods, beverages and other materials**

**Principle:**

(formate dehydrogenase)



**Kit size:** 25 assays (manual) / 250 (microplate) / 250 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 12 min

**Detection limit:** 0.0932 mg/L

**Application examples:** Wine, fruit juices, pickles, vinegar, jam, bakery products, honey, fish, meat and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by MEBAK

**Advantages**

- No wasted formate dehydrogenase solution (stable suspension supplied)
- Pyrazole incorporated to prevent alcohol dehydrogenase interference
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats




## Fructan (Hexokinase format)

Cat. No. **K-FRUCHK**

**UV-method for the determination of Fructan in foodstuffs, beverages and other materials**

**Principle:**

(sucrase + maltase)



(*exo*-inulinase + *endo*-inulinase)



(hexokinase)



(glucose-6-phosphate dehydrogenase)



(phosphoglucose isomerase)



**Kit size:** 50 assays

**Method:** Spectrophotometric at 340 nm

**Total assay time:** ~ 90 min

**Detection limit:** 1-100% of sample weight

**Application examples:** Flours, plant materials (e.g. onion), food products and other materials

**Method recognition:** This method is a modification of AOAC Method 999.03 and AACC Method 32-32.01

**Advantages**

- Very cost effective
- All reagents stable for > 12 months after preparation
- Fructan kits are available only from Megazyme
- Simple format
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standards included



Food



Feed

## Fructan (PAHBAH format)

Cat. No. **K-FRUC**

Colourimetric method for the determination of Fructan in plant products, foodstuffs and other materials

### Principle:

- (sucrase)  
**(1) Sucrose + H<sub>2</sub>O → D-glucose + D-fructose**
- (β-amylase + maltase + pullulanase)  
**(2) Starch + maltosaccharides + H<sub>2</sub>O → D-glucose**
- (borohydride)  
**(3) D-Glucose + D-fructose → D-sorbitol + D-mannitol (non-reducing)**
- (exo-inulinase + endo-inulinase)  
**(4) Fructan + H<sub>2</sub>O → D-glucose + D-fructose**
- (100°C, 6 min)  
**(5) D-Glucose + D-fructose + PAHBAH → PAHBAH colour complex**

<b>Kit size:</b>	100 assays
<b>Method:</b>	Spectrophotometric at 410 nm
<b>Total assay time:</b>	~ 90 min
<b>Detection limit:</b>	1-100% of sample weight
<b>Application examples:</b>	Flours, plant materials (e.g. onion), food products and other materials
<b>Method recognition:</b>	AOAC (Method 999.03), AACC (Method 32-32) and CODEX (Type III Method)

### Advantages

- Very cost effective
- All kit reagents stable for > 2 years after preparation
- Unaffected by high sucrose / reducing sugar concentrations
- Fructan kits are only available from Megazyme
- Simple format
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standards included



Food



Feed



Fermentation



Wine



Brewing



Dairy



Biofuels

## D-Fructose / D-Glucose

Cat. No. **K-FRUGL**

UV-method for the determination of D-Fructose and D-Glucose in foodstuffs, beverages and other materials

### Principle:

- (hexokinase)  
**(1) D-Glucose + ATP → G-6-P + ADP**
- (hexokinase)  
**(2) D-Fructose + ATP → F-6-P + ADP**
- (glucose-6-phosphate dehydrogenase)  
**(3) G-6-P + NADP<sup>+</sup> → gluconate-6-phosphate + NADPH + H<sup>+</sup>**
- (phosphoglucose isomerase)  
**(4) F-6-P ↔ G-6-P**

<b>Kit size:</b>	110 assays (manual) / 1100 (microplate) / 1100 (auto-analyser)
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Reaction time:</b>	~ 13 min
<b>Detection limit:</b>	0.66 mg/L
<b>Application examples:</b>	Wine, beer, fruit juices, soft drinks, milk, jam, honey, dietetic foods, bread, bakery products, candies, desserts, confectionery, ice-cream, fruit and vegetables, condiments, tobacco, cosmetics, pharmaceuticals, paper and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	Methods based on this principle have been accepted by AOAC, EN, NEN, NF, DIN, GOST, OIV, IFU, AIJN, MEBAK, IOCCC

### Advantages

- PVP incorporated to prevent tannin inhibition
- Validated by the University of Wine, Suze la Rousse, France
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation (manual analysis applications)
- Rapid reaction at either 25 or 37°C
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standards included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats



Food



Feed



Fermentation



Wine



Brewing



Dairy




Biofuels

## D-Fructose / D-Glucose Liquid Ready Reagents

Cat. No. **K-FRGLQR**

**UV-method suitable for auto-analyser and microplate formats for the determination of D-Fructose and D-Glucose in foodstuffs, beverages and other materials**

### Principle:



**Kit size:** 1100 assays (microplate) / 1100 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 13 min

**Detection limit:** 133 mg/L (recommended format)

**Application examples:** Wine, beer, fruit juices, soft drinks, milk, jam, honey, dietetic foods, bread, bakery products, candies, desserts, confectionery, ice-cream, fruit and vegetables, condiments, tobacco, cosmetics, pharmaceuticals, paper and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by AOAC, EN, NEN, NF, DIN, GOST, OIV, IFU, AIJN, MEBAK, IOCCC

### Advantages

- PVP incorporated to prevent tannin inhibition
- "Ready to use" liquid stable formulation
- Very competitive price (cost per test)
- All reagents stable for > 2 years
- Very rapid reaction (~ 13 min)
- Standard included
- Suitable for microplate and auto-analyser formats



Food



Feed



Fermentation



Wine



Brewing



Dairy






Biofuels

## D-Fructose / D-Glucose (MegaQuant™ format)

Cat. No. **K-FRGLMQ**

**Simple colourimetric method for the determination of D-Fructose and D-Glucose in foodstuffs, beverages and other materials**

### Principle:



**Kit size:** 60 assays

**Method:** Spectrophotometric at 505 nm

**Reaction time:** ~ 10 min

**Detection limit:** 15.4 mg/L

**Application examples:** Grape juice / must, wine, beer, fruit juices, soft drinks, milk, jam, honey, dietetic foods, bread, bakery products, candies, desserts, confectionery, ice-cream, fruit and vegetables, condiments, tobacco, cosmetics, pharmaceuticals, paper and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** *Novel method*

### Advantages

- Novel product, patented technology
- Spectrophotometer / laboratory expertise not required
- Highly stable reagents (at least three seasons use)
- Very competitive price (cost per test)
- Very simple procedure
- Rapid reaction time (~ 10 min)
- Standard included



## L-Fucose

Cat. No. **K-FUCOSE**

**UV-method for the determination of L-Fucose in plant material, polysaccharides, pharmaceuticals and other materials**

**Principle:**

(L-fucose dehydrogenase)



**Kit size:** 100 assays (manual) / 1000 (microplate) / 1020 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 10 min

**Detection limit:** 15.4 mg/L

**Application examples:** L-Fucose is present as the main component in fucoidan (a marine polysaccharide), foods, pharmaceuticals and other materials (e.g. biological samples, etc.)

**Method recognition:** **Novel method**

**Advantages**

- Very cost effective
- All reagents stable for > 2 years after preparation
- Only enzymatic kit available
- Simple format
- Rapid reaction time (~ 10 min)
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats




## Galactomannan

Cat. No. **K-GALM**

**UV-method for the determination of Galactomannan in legume seeds, foodstuffs and plant products**

**Principle:**


(β-mannanase)



(β-mannanase + α-galactosidase)



(β-galactose dehydrogenase)



**Kit size:** 100 assays

**Method:** Spectrophotometric at 340 nm

**Total assay time:** ~ 80 min

**Detection limit:** 1-100% of sample weight

**Application examples:** Seeds, milling fractions and food ingredients

**Method recognition:** **Novel method**

**Advantages**

- Galactose dehydrogenase now included in the kit
- Very cost effective
- All reagents stable for > 2 years after preparation
- Only enzymatic kit available
- Simple format
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standard included





Food



Feed



Brewing


## $\beta$ -Glucan (Mixed linkage)

Cat. No. **K-BGLU**

Colourimetric method for the determination of  $\beta$ -Glucan in cereal grains, feed, foodstuffs, beverages and other materials

### Principle:

(lichenase)



(\(\beta\)-glucosidase)



(glucose oxidase)



(peroxidase)



**Kit size:** 100 assays

**Method:** Spectrophotometric at 510 nm

**Total assay time:** ~ 100 min

**Detection limit:** 0.5-100% of sample weight

**Application examples:** Oats, barley, malt, wort, beer, food and other materials

**Method recognition:** AOAC (Method 995.16), AACC (Method 32-23.01), EBC (Methods 3.11.1, 4.16.1 and 8.11.1), ICC (Standard No. 166), RACI (Standard Method) and CODEX (Type II Method)

### Advantages

- Very cost effective
- All reagents stable for > 2 years as supplied
- Only enzymatic kit available
- Very specific
- Simple format
- Mega-Calc™ software tool is available from our website for hassle-free raw data processing
- Standards included



Food



Feed


## $\beta$ -Glucan (Yeast and mushroom)

Cat. No. **K-YBGL**

Colourimetric method for the determination of Yeast and Mushroom  $\beta$ -Glucan in yeast, mushroom, foodstuffs and other materials

### Principle:

(conc. HCl, 30°C, 45 min)



(1.3 M HCl, 100°C, 2 h)



(exo-1,3-\(\beta\)-glucanase + \(\beta\)-glucosidase)




(glucose oxidase)



(peroxidase)



(amyloglucosidase)



**Kit size:** 100 assays

**Method:** Spectrophotometric at 510 nm

**Total assay time:** ~ 100 min

**Detection limit:** 1-100% of sample weight

**Application examples:** Yeast preparations, mushroom preparations and other materials

**Method recognition:** **Novel method**

### Advantages

- Very cost effective
- All reagents stable for > 12 months after preparation
- Only enzymatic kit available
- Simple format
- Mega-Calc™ software tool is available from our website for hassle-free raw data processing
- Standards included



Food



Feed

## $\beta$ -Glucan (Yeast-Enzymatic)

Cat. No. **K-EBHLG**


Colourimetric method for the enzymatic determination of yeast  $\beta$ -Glucan and also (1-3)- $\beta$  Glucans.

### Principle:

(KOH, 4°C, 30 min)



(Gluczyme, 40°C, 16 h)



(glucose oxidase)



(peroxidase)



Kit size:	50 assays
Method:	Spectrophotometric at 510 nm
Reaction time:	~ 100 min
Detection limit:	1-100% of sample weight
Application examples:	Yeast preparations and other materials
Method recognition:	<b>Novel method</b>

### Advantages

- Very competitive price (cost per test)
- All reagents stable for > 12 months after preparation
- Mega-Calc™ software tool is available from our website for hassle-free raw data processing
- Standard included



Food



Feed



Fermentation



Brewing


## $\beta$ -Glucanase (Malt and microbial)

Cat. No. **K-MBGL**

Colourimetric method for the determination of  $\beta$ -Glucanase in malt, foodstuffs and fermentation products

### Principle:

( $\beta$ -glucanase)



Kit size:	100 assays
Method:	Spectrophotometric at 590 nm
Total assay time:	~ 30 min
Detection limit:	100 U/kg of malt
Application examples:	Malt extracts, wort, beer and other materials
Method recognition:	RACI (Standard Method)

### Advantages

- Very cost effective
- All reagents stable for > 2 years during use
- Only kit available
- Very specific
- Simple format
- Standard included



## Glucomannan

Cat. No. **K-GLUM**

**UV-method for the determination of Glucomannan in plant products, foodstuffs and other materials**

**Principle:**

( $\beta$ -mannanase)

(1) **Ac-Glucomannan + H<sub>2</sub>O → Ac-glucomanno-oligomers**

(pH 12.5)

(2) **Ac-Glucomanno-oligomers + H<sub>2</sub>O → glucomanno-oligomers + acetate**

( $\beta$ -glucosidase +  $\beta$ -mannosidase)

(3) **Glucomanno-oligomers + H<sub>2</sub>O → D-glucose + D-mannose**

(hexokinase)

(4) **D-glucose + D-mannose + ATP → G-6-P + M-6-P + ADP**

(glucose-6-phosphate dehydrogenase)

(5) **G-6-P + NADP<sup>+</sup> → gluconate-6-phosphate + NADPH + H<sup>+</sup>**

(phosphomannose isomerase)

(phosphoglucose isomerase)

(6) **M-6-P ⇌ F-6-P ⇌ G-6-P**

**Kit size:** 50 assays

**Method:** Spectrophotometric at 340 nm

**Total assay time:** 120 min

**Detection limit:** 1-100% of sample weight

**Application examples:** Jelly sweets, cosmetics, food gums and other materials

**Method recognition:** *Novel method*

### Advantages

- Very cost effective
- Only enzymatic kit available
- Simple format
- All reagents stable for > 2 years after preparation
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included



## D-Gluconic Acid / D-Glucono- $\delta$ -lactone

Cat. No. **K-GATE**

**UV-method for the determination of D-Gluconic Acid and D-Glucono- $\delta$ -lactone in foodstuffs, beverages and other materials**

**Principle:**

(gluconate kinase)

(1) **D-Gluconate + ATP → gluconate-6-phosphate + ADP**

(gluconate-6-phosphate dehydrogenase)

(2) **Gluconate-6-phosphate + NADP<sup>+</sup> → ribulose-5-phosphate + NADPH + CO<sub>2</sub> + H<sup>+</sup>**

(pH 11)

(3) **D-Glucono- $\delta$ -lactone + H<sub>2</sub>O → D-gluconate**

**Kit size:** 60 assays (manual) / 600 (microplate) / 600 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 6 min

**Detection limit:** 0.5 mg/L

**Application examples:** Wine, meat, processed meat (e.g. additives), fruit juice, dairy products, pharmaceuticals, paper and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by ISO, DIN, GOST

### Advantages

- All reagents stable for > 2 years after preparation
- Very competitive price (cost per test)
- Very rapid reaction
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats



Food



Biofuels

## D-Glucosamine (plus Glucosamine-Sulphate and -Acetate) Cat. No. K-GAMINE

UV-method for the determination of D-Glucosamine, D-Glucosamine sulphate and N-Acetyl Glucosamine in food supplements, foodstuffs, beverages and other materials

### Principle:

- (desulphation)  
**(1) Glucosamine sulphate → D-glucosamine + sulphate**  
 (deacetylation)  
**(2) N-Acetyl Glucosamine → D-glucosamine + acetate**  
 (hexokinase)  
**(3) D-Glucosamine + ATP → D-glucosamine-6-P + ATP**  
 (glucosamine 6-phosphate deaminase)  
**(4) D-Glucosamine-6-P + H<sub>2</sub>O → D-fructose-6-P + NH<sub>4</sub><sup>+</sup>**  
 (phosphoglucose isomerase)  
**(5) F-6-P ↔ G-6-P**  
 (glucose-6-phosphate dehydrogenase)  
**(6) G-6-P + NADP<sup>+</sup> → gluconate-6-phosphate + NADPH + H<sup>+</sup>**

<b>Kit size:</b>	50 assays (manual) / 500 (microplate) / 500 (auto-analyser)
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Reaction time:</b>	~ 8 min
<b>Detection limit:</b>	1.33 mg/L
<b>Application examples:</b>	Food supplements, food products and beverages
<b>Method recognition:</b>	<b>Novel method</b>



Food



Feed



Fermentation



Wine



Brewing



Dairy



Biofuels

### Advantages

- Novel product with simple format
- All reagents stable for > 2 years after preparation
- All enzymes supplied as stable suspensions
- Very rapid reaction
- Mega-Cal<sup>TM</sup> software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats

## D-Glucose (GOPOD format)

Cat. No. K-GLUC

Colourimetric method for the determination of D-Glucose in foodstuffs, beverages and other materials

### Principle:

- (glucose oxidase)  
**(1) D-Glucose + H<sub>2</sub>O + O<sub>2</sub> → D-gluconate + H<sub>2</sub>O<sub>2</sub>**  
 (peroxidase)  
**(2) 2H<sub>2</sub>O<sub>2</sub> + p-hydroxybenzoic acid + 4-aminoantipyrine → quinoneimine + 4H<sub>2</sub>O**

<b>Kit size:</b>	660 assays
<b>Method:</b>	Spectrophotometric at 510 nm
<b>Reaction time:</b>	~ 20 min
<b>Detection limit:</b>	100 mg/L
<b>Application examples:</b>	Wine, beer, fruit juices, soft drinks, milk, jam, dietetic foods, bakery products, candies, fruit and vegetables, tobacco, cosmetics, pharmaceuticals, feed, paper and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	Widely used and accepted in clinical chemistry and food analysis

### Advantages

- All reagents stable for > 12 months after preparation
- Very competitive price (cost per test)
- Simple format
- Standard included



## D-Glucose (Hexokinase format)

Cat. No. **K-GLUHK**


**UV-method for the determination of D-Glucose in foodstuffs, beverages and other materials**

**Principle:**

(hexokinase)



(glucose-6-phosphate dehydrogenase)



**Kit size:** (K-GLUHKR) 110 assays (manual) / 1100 (microplate) / 1000 (auto-analyser) or  
(K-GLUHKL) 220 assays (manual) / 2200 (microplate) / 2000 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 5 min

**Detection limit:** 0.66 mg/L

**Application examples:** Wine, beer, fruit juices, soft drinks, milk, jam, dietetic foods, bakery products, candies, fruit and vegetables, tobacco, cosmetics, pharmaceuticals (e.g. infusions), feed, paper (and cardboard) and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by AOAC, EN, NEN, NF, DIN, GOST, OIV, IFU, AIJN, MEBAK



**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Rapid reaction
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats


## Glucose Oxidase

Cat. No. **K-GLOX**

**Colourimetric method for the determination of Glucose Oxidase in foodstuffs and fermentation products**

**Principle:**

(glucose oxidase)



(peroxidase)



**Kit size:** 200 assays (manual) / 2000 (microplate) / 1960 (auto-analyser)

**Method:** Spectrophotometric at 510 nm

**Reaction time:** ~ 20 min

**Detection limit:** 10 U/L

**Application examples:** Enzyme preparations, and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** **Novel method**

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 12 months after preparation
- Simple format
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats





Feed



Biofuels


## D-Glucuronic Acid / D-Galacturonic Acid

Cat. No. **K-URONIC**


**UV-method for the determination of D-Glucuronic Acid or D-Galacturonic Acid in hydrolysates of plant material and polysaccharides and other materials**

**Principle:**

(Uronate dehydrogenase; UDH)



(Uronate dehydrogenase; UDH)



**Kit size:** 100 assays (manual) / 1000 (microplate) / 1000 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Total assay time:** ~ 10 min at 25°C or ~ 5 min at 37°C

**Detection limit:** ~ 17 mg/L

**Application examples:** Hydrolysates of plant material and polysaccharides and other materials

**Method recognition:** *Novel method*

**Advantages**

- Very cost effective
- All reagents stable for > 2 years during use
- Only test kit available
- Simple format
- *Mega-Calculator™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



Food



Feed



Fermentation



Biofuels


## α-Glucuronidase

Cat. No. **K-AGLUA**

**UV-method for the measurement of α-D-Glucuronidase in various enzyme preparations**

**Principle:**

(α-D-glucuronidase)



(uronate dehydrogenase; UDH)



**Kit size:** 50 assays (manual) / 200 (microplate)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 25 min

**Detection limit:** 17 mU/mL

**Application examples:** Enzyme preparations and other materials

**Method recognition:** *Novel method*

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 2 years as supplied
- Simple format
- *Mega-Calculator™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual and microplate formats



Food



Fermentation

## L-Glutamic Acid

Cat. No. **K-GLUT**


**Colourimetric method for the determination of L-Glutamic Acid (Monosodium Glutamate; MSG) in foodstuffs and other materials**

**Principle:**

(beef liver glutamate dehydrogenase)



(diaphorase)



**Kit size:** 60 assays (manual) / 600 (microplate) / 700 (auto-analyser)

**Method:** Spectrophotometric at 492 nm

**Reaction time:** ~ 9 min

**Detection limit:** 0.21 mg/L

**Application examples:** Fruit and vegetables (e.g. tomato), processed fruit and vegetables (e.g. tomato puree / juice, ketchup, soy sauce), condiments, processed meat products (e.g. extracts, bouillon and sausages), soup, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by ISO, GOST, NMKL

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Glutamate dehydrogenase solution stable at -20°C
- No wasted diaphorase solution (stable suspension supplied)
- Rapid reaction
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



Food



Fermentation


## L-Glutamine / Ammonia (Rapid)

Cat. No. **K-GLNAM**

**UV-method for the determination of L-Glutamine and Ammonia in cell culture media, foodstuffs and other materials**

**Principle:**

(glutaminase)



(microbial glutamate dehydrogenase)



**Kit size:** 50 assays of each

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 15 min

**Detection limit:** 0.54 mg/L (L-glutamine)  
0.06 mg/L (ammonia)

**Application examples:** Cell culture media and cultures, dietary supplements, vegetables and other materials (e.g. biological samples, etc.)

**Method recognition:** *Novel method*

**Advantages**

- Very rapid reaction due to use of high activity glutaminase and uninhibited glutamate dehydrogenase
- All enzymes supplied as stabilised suspensions
- Only enzymatic kit available
- Very cost effective
- All reagents stable for > 2 years after preparation
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included



## Glycerol

Cat. No. **K-GCROL**


**UV-method for the determination of Glycerol in foodstuffs, beverages and other materials**

**Principle:**

(glycerokinase)



(pyruvate kinase)



(L-lactate dehydrogenase)



<b>Kit size:</b>	70 assays / 700 (microplate)
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Reaction time:</b>	~ 5 min
<b>Detection limit:</b>	0.34 mg/L
<b>Application examples:</b>	Wine (and grape juice), beer, spirits, vinegar, marzipan, fruit juices, soft drinks, toothpaste, honey, tobacco, paper (and cardboard), cosmetics, pharmaceuticals, soap and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	Methods based on this principle have been accepted by OIV, MEBAK

**Advantages**

- Novel tablet format for increased stability
- Very competitive price (cost per test)
- All reagents stable for > 2 years as supplied
- Very rapid reaction
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual and microplate formats




## Glycerol (ADP-GK format)

Cat. No. **K-GCROLGK**


**UV-method for the determination of Glycerol in foodstuffs, beverages and other materials**

**Principle:**

(glycerokinase)



(ADP-GK)



(glucose-6-phosphate dehydrogenase)



<b>Kit size:</b>	70 assays (manual) / 700 (microplate) / 600 (auto-analyser)
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Reaction time:</b>	~ 7 min
<b>Detection limit:</b>	0.37 mg/L
<b>Application examples:</b>	Wine (and grape juice), beer, spirits, vinegar, marzipan, fruit juices, soft drinks, toothpaste, honey, tobacco, paper (and cardboard), cosmetics, pharmaceuticals, soap and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	<b>Novel method</b>

**Advantages**

- Novel tablet format for increased stability
- Very competitive price (cost per test)
- All reagents stable for > 2 years as supplied
- Very rapid reaction
- Positive reaction (assay proceeds with an increase in absorbance)
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



## D-3-Hydroxybutyric Acid

Cat. No. **K-HDBA**


**Colourimetric method for the determination of D-3-Hydroxybutyric Acid in foodstuffs**

**Principle:**

(3-hydroxybutyrate dehydrogenase)



(diaphorase)



<b>Kit size:</b>	60 assays (manual) / 600 (microplate) / 740 (auto-analyser)
<b>Method:</b>	Spectrophotometric at 492 nm
<b>Reaction time:</b>	~ 3 min
<b>Detection limit:</b>	0.20 mg/L
<b>Application examples:</b>	Egg, egg products (e.g. egg powder) and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	Methods based on this principle have been accepted by CEC

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Very rapid reaction (~ 3 min)
- No wasted diaphorase solution (stable suspension supplied)
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



## myo-Inositol

Cat. No. **K-INOSL**


**Colourimetric method for the determination of myo-Inositol in various sample matrices**

**Principle:**

(myo-inositol dehydrogenase)



(diaphorase)



<b>Kit size:</b>	50 assays
<b>Method:</b>	Spectrophotometric at 492 nm
<b>Reaction time:</b>	~ 10 min
<b>Detection limit:</b>	0.8 mg/L
<b>Application examples:</b>	Animal feeds, food, baby milk formulation and other materials
<b>Method recognition:</b>	<b>Novel method</b>

**Advantages**

- Very cost effective
- Reagents stable for > 2 years after preparation
- Only enzymatic kit available
- Rapid reaction
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included




## D-Isocitric Acid

Cat. No. **K-ISOC**


### UV-method for the determination of D-Isocitric Acid in foodstuffs

#### Principle:

(isocitrate dehydrogenase)



(pH 9-10)



(pH 9-10)



**Kit size:** 100 assays (manual) / 1000 (microplate) / 1000 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 3 min

**Detection limit:** 0.35 mg/L

**Application examples:** Fruit juices, fruit products, soft drinks and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by EN, NEN, NF, DIN, GOST, IFU, AIJN

#### Advantages

- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- No wasted isocitrate dehydrogenase solution (stable suspension supplied)
- Very rapid reaction
- *Mega-Cal*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats




## D-Lactic Acid

Cat. No. **K-DATE**


### UV-method for the determination of D-Lactic Acid in foodstuffs, beverages and other materials

#### Principle:

(D-lactate dehydrogenase)



(glutamate-pyruvate transaminase)



**Kit size:** 50 assays (manual) / 500 (microplate) / 450 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 5 min

**Detection limit:** 0.21 mg/L

**Application examples:** Wine, soft drinks, milk, dairy products (e.g. cream, milk / whey powder, cheese, condensed milk and yogurt), foods containing milk (e.g. dietetic foods, bakery products, baby food, chocolate, sweets and ice-cream), vinegar, fruit and vegetables, processed fruit and vegetables, meat products, food additives, paper (and cardboard), cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by DIN, GOST, IDF, EEC, EN, ISO, OIV, IFU, AIJN, MEBAK

#### Advantages

- Very rapid reaction with most samples (~ 5 min)
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- *Mega-Cal*™ software tool is available from our website for hassle-free raw data processing
- Standards included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats





Food



Fermentation



Wine



Dairy



Biofuels


## D- / L-Lactic Acid

Cat. No. **K-DLATE**


**UV-method for the determination of D-/L-Lactic Acid in foodstuffs, beverages and other materials**

**Principle:**


(D-lactate dehydrogenase)



(L-lactate dehydrogenase)



(glutamate-pyruvate transaminase)



**Kit size:**

50 assays of each

**Method:**

Spectrophotometric at 340 nm

**Reaction time:**

~ 10 min (L-lactic acid) and ~ 5 min (D-lactic acid)

**Detection limit:**

0.21 mg/L

**Application examples:**

Wine, soft drinks, milk, dairy products, foods containing milk (e.g. dietetic foods, bakery products, baby food, chocolate, sweets and ice-cream), vinegar, fruit and vegetables, processed fruit and vegetables, meat products, food additives, paper (and cardboard), cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:**

Methods based on this principle have been accepted by DIN, GOST, IDF, EEC, EN, ISO, OIV, IFU, AIJN, MEBAK

**Advantages**

- Rapid total analysis time (concurrent / flexible D- and L-lactic acid reaction format)
- D-lactate dehydrogenase reaction very rapid with most samples (~ 5 min)
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standards included
- Extended cofactors stability



Food



Fermentation



Wine



Dairy



Biofuels


## L-Lactic Acid

Cat. No. **K-LATE**


**UV-method for the determination of L-Lactic Acid in foodstuffs, beverages and other materials**

**Principle:**

(L-lactate dehydrogenase)



(glutamate-pyruvate transaminase)



**Kit size:**

50 assays (manual) / 450 (microplate) / 500 (auto-analyser)

**Method:**

Spectrophotometric at 340 nm

**Reaction time:**

~ 10 min

**Detection limit:**

0.21 mg/L

**Application examples:**

Wine, beer, soft drinks, milk, dairy products (e.g. cream, milk / whey powder, cheese, condensed milk and yogurt), foods containing milk (e.g. dietetic foods, bakery products, baby food, chocolate, sweets and ice-cream), egg, egg products (e.g. egg powder), baking additives, vinegar, fruit and vegetables, processed fruit and vegetables (e.g. tomatoes), meat products, food additives, feed, paper (and cardboard), cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:**

Methods based on this principle have been accepted by DIN, GOST, IDF, EEC, EN, ISO, OIV, IFU, AIJN, MEBAK

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Rapid reaction
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats




## Lactose / D-Galactose (Rapid)

Cat. No. **K-LACGAR**

**UV-method for the determination of Lactose and D-Galactose in foodstuffs, beverages and other materials**

**Principle:**



<b>Kit size:</b>	115 assays
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Reaction time:</b>	~ 15 min
<b>Detection limit:</b>	2.96 mg/L (lactose)
<b>Application examples:</b>	Milk, dairy products (e.g. cream, milk / whey powder, cheese, condensed milk and yogurt), foods containing milk (e.g. dietetic foods, bakery products, baby food, chocolate, sweets and ice-cream), food additives, feed, cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	Methods based on this principle have been accepted by AOAC, NBN, DIN, GOST, IDF

**Advantages**

- Very rapid reaction due to inclusion of galactose mutarotase (patented technology PCT/IE2004/00170)
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standard included




## Lactose / Sucrose / D-Glucose

Cat. No. **K-LACSU**

**Colourimetric method for the determination of Lactose, Sucrose and D-Glucose in foodstuffs, beverages and other materials**

**Principle:**



<b>Kit size:</b>	100 assays of each
<b>Method:</b>	Spectrophotometric at 510 nm
<b>Total assay time:</b>	~ 60 min
<b>Detection limit:</b>	100 mg/L
<b>Application examples:</b>	Flours, beverages, dairy products, milk, foodstuffs containing milk, cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	Used and accepted in food analysis

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 12 months after preparation
- Simple format
- Very specific
- Rapid reaction
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standards included



## Lactulose

Cat. No. **K-LACTUL**


**UV-method for the determination of Lactulose in milk and foodstuffs containing dairy products**

**Principle:**

( $\beta$ -galactosidase)




(glucose oxidase + catalase + H<sub>2</sub>O<sub>2</sub>)




(hexokinase)




(phosphoglucose isomerase)



(glucose-6-phosphate dehydrogenase)



(gluconate-6-phosphate dehydrogenase)



**Kit size:** 50 assays

**Method:** Spectrophotometric at 340 nm

**Total assay time:** ~ 120 min

**Detection limit:** 4.8 mg/L

**Application examples:** Milk, dairy products and foods containing milk

**Method recognition:** *Novel method*

### Advantages

- Twice the sensitivity of traditional hexokinase based lactulose methods
- Very cost effective
- All reagents stable for > 2 years after preparation
- *Mega-Calc*<sup>TM</sup> software tool is available from our website for hassle-free raw data processing
- Standard included




## D-Malic Acid

Cat. No. **K-DMAL**

**UV-method for the determination of D-Malic Acid in foodstuffs, beverages and other materials**

**Principle:**

(D-malate dehydrogenase)



**Kit size:** 100 assays (manual) / 1000 (microplate) / 1100 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 6 min

**Detection limit:** 0.26 mg/L

**Application examples:** Wine, beer, fruit juices, soft drinks, dietetic foods, candies, fruit and vegetables, cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by EEC, EN, DIN, OIV, IFU, and AIJN

### Advantages

- No wasted D-malate dehydrogenase solution (stable suspension supplied)
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Rapid reaction (even with difficult samples)
- *Mega-Calc*<sup>TM</sup> software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats




## L-Malic Acid

Cat. No. **K-LMAL**

**Manual format UV-method for the determination of L-Malic Acid in foodstuffs, beverages and other materials**

**Principle:**

(L-malate dehydrogenase)



(glutamate-oxaloacetate transaminase)



**Kit size:** (K-LMALR) 58 assays (manual) / 580 (microplate) or (K-LMALL) 116 assays (manual) / 1160 (microplate)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 3 min

**Detection limit:** 0.25 mg/L

**Application examples:** Wine, beer, fruit juices, soft drinks, candies, fruit and vegetables, bread, cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by AOAC, EEC, EN, NF, NEN, DIN, GOST, OIV, IFU, AIJN, MEBAK

**Advantages**

- PVP incorporated to prevent tannin inhibition
- Both enzymes supplied as stable suspensions
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Very rapid reaction (~ 3 min)
- *Mega-Cal*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual and microplate formats




## L-Malic Acid (Analyser format)

Cat. No. **K-LMALAF**

**Analyser format UV-method for the determination of L-Malic Acid in foodstuffs, beverages and other materials**

**Principle:**

(L-malate dehydrogenase)



(glutamate-oxaloacetate transaminase)



**Kit size:** 245.5 mL of prepared reagent (R1 + R2)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 3 min

**Detection limit:** 20 mg/L (recommended assay)

**Application examples:** Wine, beer, fruit juices, soft drinks, candies, fruit and vegetables, bread, cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by AOAC, EEC, EN, NF, NEN, DIN, GOST, OIV, IFU, AIJN, MEBAK

**Advantages**

- PVP incorporated to prevent tannin inhibition
- Very stable reagent when prepared for auto-analyser applications
- Linear calibration ( $R^2 \sim 0.9994$ ) up to 80 µg/mL of L-malic acid in final reaction solution
- Validated by the University of Wine, Suze la Rousse, France
- Very competitive price (cost per mL of reagent)
- Both enzymes supplied as stable suspensions
- Very rapid reaction (~ 3 min)



Food



Fermentation



Wine


## L-Malic Acid Liquid Ready Reagents

Cat. No. **K-LMALQR**


**UV-method suitable for microplate and auto-analyser formats for the determination of L-Malic Acid in foodstuffs, beverages and other materials**

**Principle:**

(L-malate dehydrogenase)



(glutamate-oxaloacetate transaminase)



**Kit size:** 1100 assays (microplate) / 1100 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 3 min

**Detection limit:** 166 mg/L (recommended format)

**Application examples:** Wine, beer, fruit juices, soft drinks, candies, fruit and vegetables, bread, cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by AOAC, EEC, EN, NF, NEN, DIN, GOST, OIV, IFU, AIJN, MEBAK

**Advantages**

- PVP incorporated to prevent tannin inhibition
- "Ready to use" liquid stable formulation
- Very competitive price (cost per test)
- All reagents stable for > 18 months
- Very rapid reaction (~ 3 min)
- Standard included
- Suitable for microplate and auto-analyser formats



Food



Fermentation



Wine


## L-Malic Acid (MegaQuant™ format)

Cat. No. **K-LMALMQ**


**Simple colourimetric method for the determination of L-Malic Acid in foodstuffs, beverages and other materials**

**Principle:**


(L-malate dehydrogenase)



(glutamate-oxaloacetate transaminase)



(diaphorase)



**Kit size:** 60 assays

**Method:** Spectrophotometric at 505 nm

**Reaction time:** ~ 6 min

**Detection limit:** 0.25 mg/L

**Application examples:** Wine, beer, fruit juices, soft drinks, candies, fruit and vegetables, bread, cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** *Novel method*

**Advantages**

- Novel product, patented technology
- Highly stable reagents (at least three seasons use)
- Very competitive price (cost per test)
- Spectrophotometer / laboratory / expertise not required
- Very simple procedure
- Rapid reaction time (~ 6 min)
- Standard included





Food



Feed



Fermentation



Brewing

## Malt Amylase

Cat. No. **K-MALTA**

Colourimetric method for the determination of  $\alpha$ -Amylase and  $\beta$ -Amylase in cereal grains, malt, food, beverages and fermentation products

### Principle:

(1)  $\alpha$ -Amylase is measured using the “Ceralpha” Method as used in K-CERA

(2)  $\beta$ -Amylase is measured using the “Betamyl-3” Method as used in K-BETA3

**Kit size:** 50 assays of each

**Method:** Spectrophotometric at 400 nm

**Total assay time:** ~ 20 min (Ceralpha Method)

~ 10 min (Betamyl-3 Method)

**Detection limit:** 0.05 U/mL

**Application examples:** Cereal flours, malts, fermentation broths and other materials

### Method recognition

**“Ceralpha” Method:** AOAC (Method 2002.01), AACC (Method 22-02.01), ICC (Standard No. 303), RACI (Standard Method) and CCFRA (Flour Testing Working Group Method 0018)

**“Betamyl-3” Method:** RACI (Standard Method)

### Advantages

- Very cost effective
- All reagents stable for > 2 years as supplied
- Only enzymatic kit available (Beta-Amylase)
- Very specific
- Simple format
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standard included



Food



Feed



Fermentation



Brewing



Biofuels

## Maltose / Sucrose / D-Glucose

Cat. No. **K-MASUG**

UV-method for the determination of Maltose, Sucrose and D-Glucose in foodstuffs, beverages and other materials

### Principle:

(1) Maltose + H<sub>2</sub>O  $\xrightarrow{(\alpha\text{-glucosidase})}$  D-glucose

(2) Sucrose + H<sub>2</sub>O  $\xrightarrow{(\beta\text{-fructosidase})}$  D-glucose + D-fructose

(3) D-Glucose + ATP  $\xrightarrow{(\text{hexokinase})}$  G-6-P + ADP

(4) G-6-P + NADP<sup>+</sup>  $\xrightarrow{(\text{glucose-6-phosphate dehydrogenase})}$  gluconate-6-phosphate + NADPH + H<sup>+</sup>

**Kit size:** 34 assays of each

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 13 min

**Detection limit:** 1.5 mg/L

**Application examples:** Beer, fruit juices, soft drinks, milk, jam, honey, dietetic foods, baby foods, bread, sugar products, bakery products, candies, desserts, confectionery, chocolate, ice-cream, fruit and vegetables, condiments, tobacco, cosmetics, pharmaceuticals, paper and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by AOAC, EN, NEN, NF, DIN, GOST, OIV, IFU, AIJN, MEBAC

### Advantages

- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Rapid reaction
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standards included




## D-Mannitol / L-Arabitol

Cat. No. **K-MANOL**


**UV-method for the determination of D-Mannitol and L-Arabitol in foodstuffs and other materials**

**Principle:**

(mannitol dehydrogenase)



(mannitol dehydrogenase)



**Kit size:** 60 assays (manual) / 600 (microplate) / 580 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 6 min

**Detection limit:** 0.50 mg/L

**Application examples:** Wine, chewing gum, dietetic foods, candies, cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** **Novel method**

**Advantages**

- Novel product (only enzymatic kit available)
- Very cost effective
- All reagents stable for > 2 years after preparation
- Simple format
- Rapid reaction
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats



## D-Mannose / D-Fructose / D-Glucose

Cat. No. **K-MANGL**


**UV-method for the determination of D-Mannose, D-Fructose and D-Glucose in foodstuffs, yeast cell preparations and other materials**

**Principle:**

(hexokinase)



(glucose-6-phosphate dehydrogenase)



(phosphomannose isomerase)

(phosphoglucose isomerase)



**Kit size:** 55 assays

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 30 min

**Detection limit:** 0.7 mg/L

**Application examples:** Foodstuffs, yeast cell preparations, enzymatic hydrolysates and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** **Novel method**

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Only enzymatic kit available
- Simple format
- Rapid reaction
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standards included



## Pectin Identification

Cat. No. **K-PECID**

**UV-method for the identification of Pectin in foodstuffs, feed and fruit juice**

**Principle:**

(pH 12.5)

**(1) Pectin + H<sub>2</sub>O → pectate + methanol**

(pectate lyase)

**(2) Pectate → 4,5-unsaturated oligogalacturonates**

**Kit size:** 500 assays

**Method:** Spectrophotometric at 235 nm

**Reaction time:** ~ 30 min

**Detection limit:** N/A

**Application examples:** Food ingredients (e.g. citrus fruit and apple) and other materials

**Method recognition:** JECFA

### Advantages

- Very cost effective
- All reagents stable for > 2 years after preparation
- Only enzymatic kit available
- Simple format
- Standards included



## Phytic Acid (Total Phosphorus)

Cat. No. **K-PHYT**

**Colourimetric method for the determination of Phytic Acid in cereal products, seed materials, animal feeds and other materials**

**Principle:**

(phytase)

**(1) Phytic acid + H<sub>2</sub>O → myo-Inositol (phosphate)<sub>n</sub> + P<sub>i</sub>**

(alkaline phosphatase)

**(2) myo-Inositol (phosphate)<sub>n</sub> → myo-inositol + P<sub>i</sub>**

**(3) P<sub>i</sub> + ammonium molybdate → 12-molybdophosphoric acid**

(diaphorase)

**(4) 12-molybdophosphoric acid + H<sub>2</sub>SO<sub>4</sub> / ascorbic acid → molybdenum blue**

**Kit size:** 50 assays

**Method:** Spectrophotometric at 655 nm

**Reaction time:** 25 min enzymic; 1 h for phosphate determination

**Detection limit:** ~ 11.3 mg phosphorus (~ 40 mg phytic acid) / 100 g material

**Application examples:** Seed materials, feeds and foodstuffs

**Method recognition:** **Novel method**

### Advantages

- Very cost effective
- All reagents stable for > 2 years after preparation
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included



## Primary Amino Nitrogen (NOPA)

Cat. No. **K-PANOPA**

**UV-method for the determination of Primary Amino Nitrogen in grape juice, must, wine and other materials**

**Principle:**

(room temperature)

**(I) Amino nitrogen + N-acetyl-L-cysteine + o-phthaldialdehyde → isoindole derivative**

**Kit size:** 100 assays (manual) / 1000 (microplate) / 1100 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 15 min

**Detection limit:** 2.59 mg N/L

**Application examples:** Grape juice, must, wine and other materials

**Method recognition:** *New method*

**Advantages**

- Simple format (absorbances read at 340 nm)
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



## Pyruvic Acid

Cat. No. **K-PYRUV**

**UV-method for the determination of Pyruvic Acid in beer, cheese, fermentation products and other materials**

**Principle:**

(D-lactate dehydrogenase)

**(I) Pyruvate + NADH + H<sup>+</sup> → D-lactic acid + NAD<sup>+</sup>**

**Kit size:** 100 assays (manual) / 1000 (microplate) / 1000 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 3 min

**Detection limit:** 0.39 mg/L

**Application examples:** Wine, beer, fruit juices, soft drinks, cheese, dietary supplements, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** *New method*

**Advantages**

- Very cost effective
- All reagents stable for > 2 years after preparation
- Very rapid reaction (~ 3 min)
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats





## Raffinose / D-Galactose

Cat. No. **K-RAFGA**

**UV-method for the determination of Raffinose (also stachyose and verbascose) and D-Galactose in legume seeds, plant materials, foodstuffs and feed**

**Principle:**



**Kit size:** 120 assays

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 40 min

**Detection limit:** 5 mg/L

**Application examples:** Cereal flours, soybean flour, by-products of sucrose manufacture and other materials

**Method recognition:** Used and accepted in food analysis

### Advantages

- Very rapid reaction due to inclusion of galactose mutarotase (patented technology)
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included





## Raffinose / Sucrose / D-Glucose

Cat. No. **K-RAFGL**

**Colourimetric method for the determination of Raffinose (also stachyose and verbascose), Sucrose and D-Glucose in legume seeds, plant materials, foodstuffs and feed**

**Principle:**



**Kit size:** 120 assays

**Method:** Spectrophotometric at 510 nm

**Reaction time:** ~ 20 min

**Detection limit:** 100 mg/L

**Application examples:** Analysis of grain legumes and other materials containing raffinose, stachyose and verbascose

**Method recognition:** Used and accepted in food analysis

### Advantages

- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Simple format
- Rapid reaction
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standards included





Feed



Biofuels


## L-Rhamnose

Cat. No. **K-RHAMNOSE**

**UV-method for the determination of L-Rhamnose in hydrolysates of plant material, polysaccharides, culture media / supernatants and other materials. Suitable for use with manual, microplate and auto-analyser formats.**

**Principle:**

(L-rhamnose dehydrogenase)



**Kit size:** 50 / 100 assays (manual) / 550 (microplate) / 550 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Total assay time:** ~ 5 min at 25°C or ~ 4 min at 37°C

**Detection limit:** ~ 1.2 mg/L

**Application examples:** Hydrolysates of plant material and polysaccharides, culture media / supernatants and other materials

**Method recognition:** *Novel method*

**Advantages**

- Very cost effective
- All reagents stable for > 2 years during use
- Only test kit available
- Simple format
- Rapid reaction (~ 5 min at 25°C or ~ 4 min at 37°C)
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



Food



Wine


## D-Sorbitol / Xylitol

Cat. No. **K-SORB**

**Colourimetric method for the determination of D-Sorbitol and Xylitol in foodstuffs and wine**

**Principle:**


(sorbitol dehydrogenase)



(sorbitol dehydrogenase)



(diaphorase)



**Kit size:** 58 assays (manual) / 580 (microplate) / 700 (auto-analyser)

**Method:** Spectrophotometric at 492 nm

**Reaction time:** ~ 15 min

**Detection limit:** 0.20 mg/L

**Application examples:** Diabetic foods (e.g. honey, jam and chocolate), dietetic foods, chewing gum, candies, fruit juice (e.g. apple juice), ice-cream, sweets, bakery products (e.g. desserts), marzipan, paper (and cardboard), cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by IFU, AIJN

**Advantages**

- Each vial of sorbitol dehydrogenase is stable for > 2 months at 4°C after dissolution
- No wasted diaphorase solution (stable suspension supplied)
- Very competitive price (cost per test)
- Reagents stable for > 2 years as supplied
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



## Starch Damage

Cat. No. **K-SDAM**

**Colourimetric method for the determination of Starch Damage in cereal flours**

**Principle:**

- (fungal  $\alpha$ -amylase)
- (1) **Damaged (or gelatinised) starch + H<sub>2</sub>O → maltodextrins**
- (amyloglucosidase)
- (2) **Maltodextrins + H<sub>2</sub>O → D-glucose**
- (glucose oxidase)
- (3) **D-Glucose + H<sub>2</sub>O + O<sub>2</sub> → D-gluconate + H<sub>2</sub>O<sub>2</sub>**
- (peroxidase)
- (4) **2H<sub>2</sub>O<sub>2</sub> + p-hydroxybenzoic acid + 4-aminoantipyrine → quinoneimine + 4H<sub>2</sub>O**

<b>Kit size:</b>	200 assays
<b>Method:</b>	Spectrophotometric at 510 nm
<b>Total assay time:</b>	~ 40 min
<b>Detection limit:</b>	0.5-100% of sample weight
<b>Application examples:</b>	Cereal flours and other materials
<b>Method recognition:</b>	AACC (Method 76-31.01), ICC (Standard No. 164), and RACI (Standard Method)

**Advantages**

- Very cost effective
- All reagents stable for > 2 years as supplied
- Only enzymatic kit available
- Very specific
- Simple format
- *Mega-Calc*<sup>TM</sup> software tool is available from our website for hassle-free raw data processing
- Standards included



## Resistant Starch

Cat. No. **K-RSTAR**

**Colourimetric method for the determination of Resistant Starch in cereal products and feeds**

**Principle:**

- ( $\alpha$ -amylase + amyloglucosidase)
- (1) **Non-resistant starch + H<sub>2</sub>O → D-glucose + maltose (trace)**
- (2) **Aqueous ethanol wash + centrifugation to remove D-glucose + maltose**
- (3) **Dissolution of resistant starch pellet in KOH and neutralisation**
- ( $\alpha$ -amylase + amyloglucosidase)
- (4) **Dissolved resistant starch + H<sub>2</sub>O → D-glucose**
- (glucose oxidase)
- (5) **D-Glucose + H<sub>2</sub>O + O<sub>2</sub> → D-gluconate + H<sub>2</sub>O<sub>2</sub>**
- (peroxidase)
- (6) **2H<sub>2</sub>O<sub>2</sub> + p-hydroxybenzoic acid + 4-aminoantipyrine → quinoneimine + 4H<sub>2</sub>O**

<b>Kit size:</b>	100 assays
<b>Method:</b>	Spectrophotometric at 510 nm
<b>Total assay time:</b>	~ 120 min (plus overnight incubation)
<b>Detection limit:</b>	2-100% of sample weight
<b>Application examples:</b>	Plant materials, starch samples and other materials
<b>Method recognition:</b>	AOAC (Method 2002.02), AACC (Method 32-40.01) and CODEX (Type II Method)

**Advantages**

- Very cost effective
- All reagents stable for > 2 years as supplied
- Only enzymatic kit available
- Measures enzyme resistant starch
- Simple format
- *Mega-Calc*<sup>TM</sup> software tool is available from our website for hassle-free raw data processing
- Standard included



Food



Feed



Fermentation



Wine



Brewing

## Total Starch (GOPOD format)

Cat. No. K-TSTA

**Colourimetric method for the determination of Total Starch in cereal products, feeds, foodstuffs and other materials**

**Principle:**

( $\alpha$ -amylase, 100°C  $\pm$  DMSO)

(1) **Starch granules + H<sub>2</sub>O  $\rightarrow$  maltodextrins**

(amyloglucosidase)

(2) **Maltodextrins + H<sub>2</sub>O  $\rightarrow$  D-glucose**

(glucose oxidase)

(3) **D-Glucose + H<sub>2</sub>O + O<sub>2</sub>  $\rightarrow$  D-gluconate + H<sub>2</sub>O<sub>2</sub>**

(peroxidase)

(4) **2H<sub>2</sub>O<sub>2</sub> + *p*-hydroxybenzoic acid + 4-aminoantipyrine  $\rightarrow$  quinoneimine + 4H<sub>2</sub>O**

**Kit size:** 100 assays

**Method:** Spectrophotometric at 510 nm

**Total assay time:** ~ 90 min

**Detection limit:** 1-100% of sample weight

**Application examples:** Cereal flours, food products and other materials

**Method recognition:** AOAC (Method 996.11), AACC (Method 76-13.01), ICC (Standard Method No. 168), and RACI (Standard Method)

### Advantages

- Very competitive price (cost per test)
- All reagents stable for > 12 months after preparation
- Simple format
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included



Food



Feed



Fermentation



Wine



Brewing

## Total Starch (Hexokinase format)

Cat. No. K-TSHK

**UV-method for the determination of Total Starch in grains, animal feeds, foodstuffs and other materials**

**Principle:**

( $\alpha$ -amylase, 100°C  $\pm$  DMSO)

(1) **Starch granules + H<sub>2</sub>O  $\rightarrow$  maltodextrins**

(amyloglucosidase)

(2) **Maltodextrins + H<sub>2</sub>O  $\rightarrow$  D-glucose**

(hexokinase)

(3) **D-Glucose + ATP  $\rightarrow$  G-6-P + ADP**

(glucose-6-phosphate dehydrogenase)

(4) **G-6-P + NADP<sup>+</sup>  $\rightarrow$  gluconate-6-phosphate + NADPH + H<sup>+</sup>**

**Kit size:** 100 assays

**Method:** Spectrophotometric at 340 nm

**Total assay time:** ~ 90 min

**Detection limit:** 1-100% of sample weight

**Application examples:** Cereal flours, food products and other materials

**Method recognition:** AOAC (Method 996.11), AACC (Method 76-13.01), ICC (Standard Method No. 168), and RACI (Standard Method)

### Advantages

- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Simple format
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included



## Succinic Acid

Cat. No. **K-SUCC**


**UV-method for the determination of Succinic Acid in foodstuffs, feed, wine and other materials**

**Principle:**

(succinyl-CoA synthetase)



(pyruvate kinase)



(L-lactate dehydrogenase)



**Kit size:** 20 assays (manual) / 200 (microplate) / 270 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 6 min

**Detection limit:** 0.26 mg/L

**Application examples:** Wine, fruit and vegetables, soy sauce, cheese, egg, egg products and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by EEC

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 2 years as supplied
- Very rapid reaction (even at room temperature)
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats



## Sucrose / D-Fructose / D-Glucose

Cat. No. **K-SUFRG**

**UV-method for the determination of Sucrose, D-Fructose and D-Glucose in foodstuffs, beverages and other materials**

**Principle:**

(β-fructosidase)




(hexokinase)




(hexokinase)



(glucose-6-phosphate dehydrogenase)



(phosphoglucose isomerase)



**Kit size:** 50 assays of each

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 30 min

**Detection limit:** 1.38 mg/L

**Application examples:** Beer, fruit juices, soft drinks, milk, jam, honey, dietetic foods, bread, bakery products, dairy products, candies, desserts, confectionery, sweets, ice-cream, fruit and vegetables (e.g. potato), meat products (e.g. sausage), condiments (e.g. ketchup and mustard), feed, tobacco, cosmetics, pharmaceuticals, paper and other materials

**Method recognition:** Methods based on this principle have been accepted by NF, EN, NEN, DIN, GOST, IFU, AIJN, MEBAK, IOCCC

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- Rapid reaction
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Stabilised D-glucose / D-fructose standard solution included
- Extended cofactors stability



Food



Feed



Fermentation



Wine



Biofuels

## Sucrose / D-Glucose (GOPOD format)

Cat. No. K-SUCGL

**Colourimetric method for the determination of Sucrose and D-Glucose in foodstuffs, beverages and other materials**

**Principle:**

(glucose oxidase)



(peroxidase)



( $\beta$ -fructosidase)



**Kit size:** 250 assays

**Method:** Spectrophotometric at 510 nm

**Reaction time:** ~ 30 min

**Detection limit:** 100 mg/L

**Application examples:** Beer, fruit juices, soft drinks, coffee, milk, jam, honey, dietetic foods, bread, bakery products, candies, chocolate, desserts, confectionery, ice-cream, fruit and vegetables, condiments, tobacco, cosmetics, pharmaceuticals, paper and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Used and accepted in food analysis

**Advantages**

- Very competitive price (cost per test)
- All reagents stable for > 12 months after preparation
- Simple format
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standards included



Food



Wine

## Total Sulphite

Cat. No. K-TSULPH

**Colourimetric methods for the determination of Total Sulphite in wine, fruit juice, foodstuffs and other materials**

**Principle:**

The Total Sulphite assay is based on the reaction principle between thiol groups and Ellman's reagent

**Kit size:** 80 assays (manual) / 800 (microplate) / 800 (auto-analyser)

**Method:** Spectrophotometric at 405 nm

**Total assay time:** ~ 6 min

**Detection limit:** ~ 5 mg/L

**Application examples:** Wine, fruit juice, sea food, food stuffs and other materials

**Method recognition:** Validated for red and white wines at the Bundesamt für Weinbau, Austria.  
Used widely in the wine industry

**Advantages**

- "Ready to use" liquid stable formulation
- Very competitive price (cost per test)
- All reagents stable for > 18 months
- Very rapid reaction
- *Mega-Calc*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats






## Total Sulphite (Enzymatic)

Cat. No. **K-ETSULPH**

**UV-method for the determination of Total Sulphite (SO<sub>2</sub><sup>2-</sup>) in beverages, foodstuffs and other materials**

### Principle:



<b>Kit size:</b>	50 assays (manual) / 500 (microplate) / 588 (auto-analyser)
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Total assay time:</b>	~ 30 min
<b>Detection limit:</b>	0.3 mg/L
<b>Application examples:</b>	Wine, beer, fruit juices, soft drinks, jam, fruit and vegetables, pharmaceuticals and other materials
<b>Method recognition:</b>	Methods based on this principle have been accepted by DIN, EN, MEBAK and NMKL

### Advantages

- Very cost effective
- All reagents stable for > 2 years during use
- Simple format
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability
- Suitable for manual, microplate and auto-analyser formats



## Total and Free Sulphite

Cat. No. **K-SULPH**

**Colourimetric methods for the determination of Total and Free Sulphite in wine, fruit juice, foodstuffs and other materials**

### Principle:

**The Total Sulphite assay is based on the reaction principle between thiol groups and Ellman's reagent**

**The Free Sulphite assay is based on the reaction principle of SO<sub>2</sub>, fuchsin and aldehyde binding**

<b>Kit size:</b>	40 assays (manual) / 400 (microplate) / 400 (auto-analyser)
<b>Method:</b>	Total sulphite: Spectrophotometric at 405 nm Free sulphite: Spectrophotometric at 575 nm
<b>Total assay time:</b>	Total sulphite: ~ 6 min Free sulphite: ~ 9 min
<b>Detection limit:</b>	Total sulphite: ~ 5 mg/L Free sulphite: ~ 2 mg/L
<b>Application examples:</b>	Wine, fruit juice, seafood, food stuffs and other materials
<b>Method recognition:</b>	Validated for red and white wines at the Bundesamt für Weinbau, Austria. Used widely in the wine industry

### Advantages

- "Ready to use" liquid stable formulation
- Very competitive price (cost per test)
- All reagents stable for > 18 months
- Very rapid reaction
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



## Tartaric Acid

Cat. No. **K-TART**

**Colourimetric method for the determination of Tartaric Acid in wine, fruit juice and other materials**

### Principle:

The Tartaric acid assay is based on the reaction principles between tartaric acid and vanadate

<b>Kit size:</b>	200 assays (manual) / 2000 (microplate) / 2000 (auto-analyser)
<b>Method:</b>	Spectrophotometric at 505 nm
<b>Total assay time:</b>	~ 4 min
<b>Detection limit:</b>	~ 108 mg/L
<b>Application examples:</b>	Wine, fruit juice and other materials
<b>Method recognition:</b>	Used widely in the wine industry

### Advantages

- “Ready to use” liquid stable formulation
- Very competitive price (cost per test)
- All reagents stable for > 1 year
- Very rapid reaction
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats





## Trehalose / D-Glucose

Cat. No. **K-TREH**

**UV-method for the determination of Trehalose and D-Glucose in foodstuffs, beverages, and other materials**

### Principle:



<b>Kit size:</b>	100 assays (manual) / 1000 (microplate) / 1100 (auto-analyser)
<b>Method:</b>	Spectrophotometric at 340 nm
<b>Reaction time:</b>	~ 8 min
<b>Detection limit:</b>	37.5 mg/L
<b>Application examples:</b>	Honey, mushrooms, bread, beer, seafood (e.g. lobster and shrimp), fruit juices, purees and fillings, nutrition bars, surimi, dehydrated fruits and vegetables, fruit products, white chocolate, sports drinks, dairy products, egg products, soups and sauces, confectionery, chewing gum, cosmetics, pharmaceuticals and other materials (e.g. biological cultures, samples, etc.)
<b>Method recognition:</b>	<b>Novel method</b>

### Advantages

- Only enzymatic kit available
- Very cost effective
- All reagents stable for > 2 years after preparation
- Very rapid reaction
- *Mega-Calc™* software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



Food



Feed



Fermentation



Wine





Dairy

## Urea / Ammonia (Rapid)

Cat. No. **K-URAMR**

**UV-method for the determination of Urea and Ammonia in foodstuffs, beverages and other materials**

**Principle:**



**Kit size:** 50 assays of each

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 8 min

**Detection limit:** 0.13 mg/L (urea)  
0.07 mg/L (ammonia)

**Application examples:** Wine, grape juice, must, fruit juices, soft drinks, milk, cheese, meat, processed meat, bakery products, seafood, fertilizers, feed, pharmaceuticals, cosmetics, water (e.g. swimming-pool water), Kjeldahl analysis, paper (and cardboard) and other materials (e.g. biological cultures, samples, etc.)

**Method recognition:** Methods based on this principle have been accepted by NEN, MEBAK

**Advantages**

- Very rapid reaction due to use of uninhibited glutamate dehydrogenase
- Enzymes supplied as stable suspensions
- Very competitive price (cost per test)
- All reagents stable for > 2 years after preparation
- *Mega-Calculator*™ software tool is available from our website for hassle-free raw data processing
- Standard included
- Extended cofactors stability



Feed




Biofuels

## Xylanase (Azo-Wax format)

Cat. No. **K-AZOWAX**

**Colourimetric method for the determination of Xylanase in feed, foodstuffs and other materials**

**Principle:**



**Kit size:** 200 assays

**Method:** Based on use of Azo-WAX reagent (590 nm)

**Total assay time:** ~ 45 min

**Detection limit:** 0.2 U/mL of assay solution

**Application examples:** Animal feeds, enzyme preparations, bread improver mixtures and other materials

**Method recognition:** Used widely in the feed industry

**Advantages**

- Very cost effective
- All reagents stable for > 2 years
- Only test kit available
- Simple format
- Standard included



Feed



Biofuels

## Xylanase (Xylazyme AX format)

Cat. No. **K-XYLS**

**Colourimetric method for the determination of Xylanase in feed, foodstuffs and other materials**

**Principle:**

( $\beta$ -xylanase)

**(1) Xylazyme AX (water insoluble) + H<sub>2</sub>O → water soluble dyed xylan fragments**

**Kit size:** 200 assays

**Method:** Based on use of Xylazyme AX tablets (590 nm)

**Total assay time:** ~ 45 min

**Detection limit:** 0.02 U/mL of assay solution

**Application examples:** Animal feeds, enzyme preparations, bread improver mixtures and other materials

**Method recognition:** Used widely in the feed industry

**Advantages**

- Very cost effective
- All reagents stable for > 2 years during use
- Only test kit available
- Simple format
- Standards included



Food



Feed



Fermentation



Biofuels

## D-Xylose

Cat. No. **K-XYLOSE**

**UV-method for the determination of D-Xylose in fermentation broths and hydrolysates of plant material and polysaccharides**

**Principle:**

(xylose mutarotase)

**(1)  $\alpha$ -D-Xylose ↔  $\beta$ -D-xylose**

( $\beta$ -xylose dehydrogenase)

**(2)  $\beta$ -D-Xylose + NAD<sup>+</sup> → D-xylonic acid + NADH + H<sup>+</sup>**

**Kit size:** 100 assays (manual) / 1000 (microplate) / 1300 (auto-analyser)

**Method:** Spectrophotometric at 340 nm

**Reaction time:** ~ 6 min

**Detection limit:** 0.7 mg/L

**Application examples:** Analysis of D-xylose in fermentation broths and hydrolysates of plant material and polysaccharides

**Method recognition:** **Novel method**

**Advantages**

- Very cost effective
- Reagents stable for > 2 years after preparation
- Only enzymatic kit available
- Rapid reaction (~ 6 min)
- *Mega-Calc*<sup>TM</sup> software tool is available from our website for hassle-free raw data processing
- Standard included
- Suitable for manual, microplate and auto-analyser formats



Food



Feed



Fermentation



Brewing


## Limit-Dextrinase / Pullulanase (Under Development)

Cat. No. K-LDPU

Available soon...

Colourimetric method for the determination of Limit-Dextrinase or Pullulanase in grain extracts, feed and fermentation products

### Principle:



**Kit size:** 100 / 200 assays  
**Method:** Spectrophotometric at 400 nm (BzCNP $G_3G_3$ )  
 Fluorimetric at 450nm (BzMUG $G_3G_3$ )  
**Total assay time:** ~ 10 min  
**Detection limit:** 0.05 u/mL  
**Application examples:** Cereal flours, fermentation broths, and other materials  
**Method recognition:** *Novel method*

### Advantages

- Suitable for manual and auto-analyser formats
- Very cost effective
- All reagents stable for > 2 years after preparation
- Very specific
- Simple format
- Mega-Calc™ software tool is available from our website for hassle-free raw data processing
- Standard included



A close-up photograph of a person wearing white gloves holding a glass flask. The neck of the flask is wrapped in crinkled aluminum foil and contains a white, fluffy substance. The main body of the flask is filled with a translucent orange liquid. The background is blurred, showing laboratory equipment.

ENZYMES...

# Enzymes

Megazyme stocks a wide range of ultra-pure enzymes for use in analytical, diagnostic and research applications. These enzymes are either produced recombinantly through various expression systems or are purified using conventional protein purification techniques from crude industrial enzyme formulations.

## Diversity

Megazyme's enzyme portfolio is constantly expanding and we are always keen to receive enzyme suggestions from customers. The bulk of Megazyme's enzymes are carbohydrate acting enzymes (CAZymes) and cover 36 CAZy families - 31 glycoside hydrolase (GH) families, three polysaccharide lyase (PL) families and two carbohydrate esterase (CE) families. A wide range of analytical enzymes including kinases and isomerases are also available.

## Purity


Megazyme prides itself on the purity of our enzymes and reagents. As an example we supply the enzymes required for the analysis of dietary fiber by the official AOAC methods, namely  $\alpha$ -amylase, amyloglucosidase and protease. The quality of these enzymes is of paramount importance for use in this assay and Megazyme proudly provides the highest purity level commercially available globally.

## Stability


Megazyme enzymes are provided in various formats including freeze-dried powder, ammonium sulphate suspension and 50% glycerol solution. All enzymes have undergone rigorous stability studies and recommended storage conditions are provided with every product.

## Background Information


Every enzyme is extensively characterised by our R&D scientists. A data sheet is shipped with each enzyme describing specific activity on a suitable substrate and relative activity on any other relevant substrates. The data sheet also describes pH and temperature activity optima along with pH and temperature stability data. These data sheets are also available on every product page at [www.megazyme.com](http://www.megazyme.com).



$\alpha$ -Amylase



$\alpha$ -Galactosidase



Invertase

Cat. No.	Product
<b>ENZYMES</b>	
E-ACSBS	<sup>Rec</sup> Acetyl-CoA synthetase ( <i>B. subtilis</i> )
E-OGLEYF	<sup>Rec</sup> endo- $\alpha$ -N-Acetylgalactosaminidase ( <i>E. faecalis</i> )
E-ANAGM	<sup>Rec</sup> $\alpha$ -N-Acetylgalactosaminidase (microbial)
E-BNAHEX	<sup>Rec</sup> $\beta$ -N-Acetylhexosaminidase (microbial)
E-AXEAO	<sup>Rec</sup> Acetylxyylan esterase ( <i>Orpinomyces</i> sp.)
E-ACPEC	<sup>Rec</sup> Phosphatase (acid) ( <i>E. coli</i> )
E-AMPK	<sup>Rec</sup> Adenylate kinase (myokinase) (prokaryote)
E-ADHEC	<sup>Rec</sup> Alcohol dehydrogenase ( <i>E. coli</i> )
E-ALGLS	<sup>Rec</sup> Alginate lyase ( <i>Sphingomonas</i> sp.)
E-ALPEC	<sup>Rec</sup> Phosphatase (Alkaline) ( <i>E. coli</i> )
E-ANAAM	$\alpha$ -Amylase ( <i>A. oryzae</i> )
E-BAASS	$\alpha$ -Amylase ( <i>B. amyloliquefaciens</i> )
E-BLAAM	$\alpha$ -Amylase ( <i>B. licheniformis</i> )
E-BAM	$\alpha$ -Amylase ( <i>B. licheniformis</i> ) <sup>NEW</sup>
E-PANAA	$\alpha$ -Amylase (Porcine Pancreatic)
E-BARBL	$\beta$ -Amylase (barley; liquid)
E-BARBP	$\beta$ -Amylase (barley; powder)
E-BAMBC	<sup>Rec</sup> $\beta$ -Amylase ( <i>B. cereus</i> )
E-AMGDF	Amyloglucosidase ( <i>A. niger</i> ) <sup>NEW</sup>

Cat. No.	Product
E-AMGDFPD	Amyloglucosidase (powder) ( <i>A. niger</i> ) <sup>NEW</sup>
E-AMGFR	Amyloglucosidase ( <i>A. niger</i> )
E-AMGPU	Amyloglucosidase ( <i>Rhizopus</i> sp.)
E-EARAB	endo-1,5- $\alpha$ -L-Arabinanase ( <i>A. niger</i> )
E-AFASE	$\alpha$ -L-Arabinofuranosidase ( <i>A. niger</i> )
E-AFAM2	<sup>Rec</sup> $\alpha$ -L-Arabinofuranosidase (novel specificity)
E-ABFCJ	<sup>Rec</sup> $\alpha$ -L-Arabinofuranosidase ( <i>C. japonicus</i> )
E-ABFCT	<sup>Rec</sup> $\alpha$ -L-Arabinofuranosidase ( <i>C. thermocellum</i> )
E-ABFUM	<sup>Rec</sup> $\alpha$ -L-Arabinofuranosidase ( <i>U. maydis</i> ) <sup>NEW</sup>
E-ARBACJ	<sup>Rec</sup> endo-/exo-Arabinanase ( <i>C. japonicus</i> )
E-ASNEC	<sup>Rec</sup> Asparaginase ( <i>E. coli</i> )
E-DIPEP	<sup>Rec</sup> $\alpha$ -Aspartyl dipeptidase ( <i>E. coli</i> )
E-CBHI	Cellobiohydrolase I ( <i>T. longibrachiatum</i> )
E-CELAN	Cellulase (endo-1,4- $\beta$ -D-glucanase) ( <i>A. niger</i> )
E-CELBA	<sup>Rec</sup> Cellulase (endo-1,4- $\beta$ -D-glucanase) ( <i>B. amyloliquefaciens</i> )
E-CELTE	Cellulase (endo-1,4- $\beta$ -D-glucanase) ( <i>T. emersonii</i> )
E-CELTR	Cellulase (endo-1,4- $\beta$ -D-glucanase) ( <i>T. longibrachiatum</i> )
E-CELTMT	<sup>Rec</sup> Cellulase (endo-1,4- $\beta$ -D-glucanase) ( <i>T. maritima</i> )
E-CITEC	<sup>Rec</sup> Citrate synthase ( <i>E. coli</i> )
E-CREA	<sup>Rec</sup> Creatinase ( <i>Bacillus</i> sp.)



Cat. No.	Product
E-CMPK	<sup>Rec</sup> Cytidylate kinase (prokaryote)
E-DIAEC	<sup>Rec</sup> Diaphorase ( <i>E. coli</i> )
E-FAERU	<sup>Rec</sup> Feruloyl esterase (rumen microorganism)
E-FAEZCT	<sup>Rec</sup> Feruloyl esterase ( <i>C. thermocellum</i> )
E-FDHCB	<sup>Rec</sup> Formate dehydrogenase ( <i>C. boidinii</i> )
E-FRMXLQ	Fructanase mixture (purified; liquid)
E-FRMXPD	Fructanase mixture (purified; powder)
E-FUCTM	<sup>Rec</sup> $\alpha$ -Fucosidase (thermostable) ( <i>T. maritima</i> ) <sup>NEW</sup>
E-FUCM	1,2- $\alpha$ -L-Fucosidase (microbial)
E-EGALN	<i>endo</i> -1,4- $\beta$ -D-Galactanase ( <i>A. niger</i> )
E-GALCJ	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -D-Galactanase ( <i>C. japonicus</i> )
E-GALCT	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -D-Galactanase ( <i>C. thermocellum</i> )
E-GALDH	<sup>Rec</sup> Galactose dehydrogenase (soil prokaryote)
E-GALMUT	<sup>Rec</sup> Galactose dehydrogenase / Galactose mutarotase
E-AGLAN	$\alpha$ -Galactosidase ( <i>A. niger</i> )
E-AGLANP	$\alpha$ -Galactosidase ( <i>A. niger</i> ) powder
E-AGALPS	<sup>Rec</sup> $\alpha$ -Galactosidase ( <i>P. simplicissimum</i> ) <sup>NEW</sup>
E-AGLGU	$\alpha$ -Galactosidase (guar)
E-BGLAN	$\beta$ -Galactosidase ( <i>A. niger</i> )
E-LICACT	<sup>Rec</sup> Non-specific <i>endo</i> -1,3(4)- $\beta$ -Glucanase ( <i>C. thermocellum</i> )
E-EXBGOS	<i>exo</i> -1,3- $\beta$ -D-Glucanase / $\beta$ -Glucosidase
E-LAMSE	<i>endo</i> -1,3- $\beta$ -D-Glucanase ( <i>Trichoderma</i> sp.)
E-LAMHV	<sup>Rec</sup> <i>endo</i> -1,3- $\beta$ -D-Glucanase (barley)
E-EXG5AO	<sup>Rec</sup> <i>exo</i> -1,3- $\beta$ -D-Glucanase ( <i>Aspergillus oryzae</i> )
E-EXBGL	<i>exo</i> -1,3- $\beta$ -D-Glucanase ( <i>Trichoderma</i> sp.)
E-EXBGTV	<sup>Rec</sup> <i>exo</i> -1,3- $\beta$ -D-Glucanase ( <i>T. virens</i> )
E-GAMP	<sup>Rec</sup> Glucoamylase P ( <i>H. resinae</i> ) <sup>NEW</sup>
E-GLUKEC	<sup>Rec</sup> Glucokinase ( <i>E. coli</i> )
E-GOX	Glucose oxidase ( <i>Aspergillus</i> sp.) <sup>NEW</sup>
E-GOXCA	Glucose oxidase / Catalase mixture (eukaryote)
E-GPDH5	Glucose-6-phosphate dehydrogenase ( <i>L. mesenteroides</i> )
E-GPDHEC	<sup>Rec</sup> Glucose-6-phosphate dehydrogenase ( <i>E. coli</i> )
E-TSAGL	$\alpha$ -Glucosidase ( <i>B. stearothermophilus</i> )
E-TSAGS	<sup>Rec</sup> $\alpha$ -Glucosidase ( <i>B. stearothermophilus</i> )
E-AGLUTM	<sup>Rec</sup> $\alpha$ -Glucosidase (thermostable) ( <i>T. maritima</i> ) <sup>NEW</sup>
E-MALTS	$\alpha$ -Glucosidase (maltase) (yeast)
E-TRNGL	$\alpha$ -Glucosidase (transglucosidase) ( <i>A. niger</i> )
E-OAGUM	Oligo- $\alpha$ -1,6-Glucosidase
E-MALBS	Oligo- $\alpha$ -1,4(6)-Glucosidase ( <i>B. subtilis</i> ) <sup>NEW</sup>
E-BGLUC	$\beta$ -Glucosidase ( <i>A. niger</i> )
E-BGOSAG	<sup>Rec</sup> $\beta$ -Glucosidase ( <i>Agrobacterium</i> sp.)
E-BGOSPC	<sup>Rec</sup> $\beta$ -Glucosidase (thermostable) ( <i>P. chrysosporium</i> )
E-BGOSTM	<sup>Rec</sup> $\beta$ -Glucosidase (thermostable) ( <i>T. maritima</i> )
E-BGLAEC	<sup>Rec</sup> $\beta$ -Glucuronidase ( <i>E. coli</i> )
E-AGUBS	<sup>Rec</sup> $\alpha$ -Glucuronidase ( <i>G. stearothermophilus</i> )
E-GOTEC	<sup>Rec</sup> Glutamate oxaloacetate transaminase ( <i>E. coli</i> )
E-GPTBS	<sup>Rec</sup> Glutamate pyruvate transaminase ( <i>B. subtilis</i> )
E-GLUTEC	<sup>Rec</sup> Glutaminase ( <i>E. coli</i> )
E-GPO	<sup>Rec</sup> Glycerol 3-phosphate oxidase
E-GMPK	<sup>Rec</sup> Guanylate kinase (prokaryote)
E-HEX10	Hexokinase (yeast)
E-HKGDH	Hexokinase / Glucose-6-phosphate dehydrogenase
E-HYLSP	<sup>Rec</sup> Hyaluronate lyase (novel specificity) (soil prokaryote)
E-HBDH	<sup>Rec</sup> 3-Hydroxybutyrate dehydrogenase (prokaryote)
E-INDHBS	<sup>Rec</sup> <i>myo</i> -Inositol dehydrogenase ( <i>B. subtilis</i> )
E-ENDOIAN	<sup>Rec</sup> <i>endo</i> -Inulinase ( <i>A. niger</i> )
E-EXOIAN	<sup>Rec</sup> <i>exo</i> -Inulinase ( <i>A. niger</i> )
E-INVRT	Invertase (fructofuranosidase) (yeast)

Cat. No.	Product
E-INVDP	Invertase (powder)
E-ISAMY	Isoamylase (glycogen 6-glucanohydrolase)
E-ICDHBS	<sup>Rec</sup> Isocitrate dehydrogenase ( <i>B. subtilis</i> )
E-DLDHLM	<sup>Rec</sup> D-Lactate dehydrogenase ( <i>L. mesenteroides</i> )
E-LLDHP	L-Lactate dehydrogenase (Porcine)
E-LICHN	Lichenase ( <i>endo</i> -1,3(4)- $\beta$ -D-glucanase) ( <i>Bacillus</i> sp.)
E-DMDHEC	<sup>Rec</sup> D-Malate dehydrogenase ( <i>E. coli</i> )
E-LMDHEC	<sup>Rec</sup> L-Malate dehydrogenase ( <i>E. coli</i> )
E-BMANN	<i>endo</i> -1,4- $\beta$ -Mannanase ( <i>A. niger</i> )
E-BMABS	<i>endo</i> -1,4- $\beta$ -Mannanase ( <i>Bacillus</i> sp.)
E-BMACJ	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -Mannanase ( <i>C. japonicus</i> )
E-BMATM	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -Mannanase ( <i>T. maritima</i> )
E-BMABC	<sup>Rec</sup> $\beta$ -Mannanase ( <i>B. circulans</i> ) <sup>NEW</sup>
E-MNHHPF	<sup>Rec</sup> Mannitol dehydrogenase ( <i>P. fluorescens</i> )
E-BMOSCF	<sup>Rec</sup> $\beta$ -Mannosidase ( <i>C. fimi</i> )
E-MAST	Malt Amylase Standard
E-PEROX	Peroxidase (Horse radish) <sup>NEW</sup>
E-NADHPO	<sup>Rec</sup> Peroxidase (NADH peroxidase) ( <i>E. faecalis</i> )
E-PCLYAN	Pectate lyase ( <i>Aspergillus</i> sp.)
E-PCLYAN2	Pectate lyase ( <i>Aspergillus</i> sp.)
E-PLYCJ	<sup>Rec</sup> Pectate lyase ( <i>C. japonicus</i> )
E-PGDHEC	<sup>Rec</sup> 6-Phosphogluconate dehydrogenase ( <i>E. coli</i> )
E-PGIBS	<sup>Rec</sup> Phosphoglucose isomerase ( <i>B. subtilis</i> )
E-PGIEC	<sup>Rec</sup> Phosphoglucose isomerase ( <i>E. coli</i> )
E-PGISC	<sup>Rec</sup> Phosphoglucose isomerase ( <i>S. cerevisiae</i> )
E-PGM	<sup>Rec</sup> Mutase ( $\alpha$ -Phosphoglucomutase) (microbial)
E-PMIEC	<sup>Rec</sup> Phosphomannose isomerase ( <i>E. coli</i> )
E-PTABS	<sup>Rec</sup> Phosphotransacetylase ( <i>B. subtilis</i> )
E-PGALUSP	<i>endo</i> -Polygalacturonanase M2 ( <i>A. aculeatus</i> )
E-BSPRPD	Protease (subtilisin A) ( <i>B. licheniformis</i> )
E-BSPRT	Protease (subtilisin A) ( <i>B. licheniformis</i> )
E-PULKP	Pullulanase M1 ( <i>K. planticola</i> )
E-PULBL	Pullulanase M2 ( <i>B. licheniformis</i> )
E-ISPUAN	Isopullulanase ( <i>A. niger</i> ) <sup>NEW</sup>
E-RHAMS	<sup>Rec</sup> $\alpha$ -Rhamnosidase (prokaryote)
E-SIALCP	<sup>Rec</sup> <i>exo</i> - $\alpha$ -Sialidase ( <i>C. perfringens</i> )
E-SIALST	<sup>Rec</sup> <i>exo</i> - $\alpha$ -Sialidase ( <i>S. typhimurium</i> )
E-SCOAS	<sup>Rec</sup> Succinyl-CoA synthetase (prokaryote)
E-SUCR	Sucrase (maltase) (yeast)
E-SUCRBG	Sucrase plus $\beta$ -Galactosidase
E-TMPK	<sup>Rec</sup> Thymidylate kinase (prokaryote)
E-TREH	<sup>Rec</sup> Trehalase (prokaryote)
E-UAO	<sup>Rec</sup> Uricase (eukaryote)
E-XANLB	<sup>Rec</sup> Xanthan lyase ( <i>Bacillus</i> sp.)
E-XYTRI	<i>endo</i> -1,4- $\beta$ -Xylanase M1 ( <i>T. viride</i> )
E-XYTR3	<i>endo</i> -1,4- $\beta$ -Xylanase M3 ( <i>T. longibrachiatum</i> ; pl 9.0)
E-XYAN4	<i>endo</i> -1,4- $\beta$ -Xylanase M4 ( <i>A. niger</i> )
E-XYRU6	<i>endo</i> -1,4- $\beta$ -Xylanase M6 (rumen microorganism)
E-XYNAP	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -Xylanase ( <i>A. punctata</i> )
E-XYNBS	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -Xylanase ( <i>B. stearothermophilus</i> T6)
E-XYNACJ	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -Xylanase ( <i>C. japonicus</i> )
E-XYNBCM	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -Xylanase ( <i>C. mixtus</i> )
E-XYLNP	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -Xylanase ( <i>N. patriciarum</i> )
E-XGP74	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -Glucanase ( <i>Paenibacillus</i> sp.)
E-XYLATM	<sup>Rec</sup> <i>endo</i> -1,4- $\beta$ -Xylanase ( <i>T. maritima</i> )
E-XEGP	<sup>Rec</sup> Xyloglucanase ( <i>Paenibacillus</i> sp.)
E-BXSEBP	<sup>Rec</sup> <i>exo</i> -1,4- $\beta$ -Xylosidase ( <i>B. pumilus</i> )
E-BXSR	<sup>Rec</sup> <i>exo</i> -1,4- $\beta$ -Xylosidase ( <i>S. ruminantium</i> ; regular)
<sup>Rec</sup>	(recombinant enzyme)

A scientist wearing a white lab coat, safety glasses, and blue gloves is working in a laboratory. The scientist is holding a small vial and adding its contents to a flask on a magnetic stirrer. The flask contains a blue liquid. In the background, there are various laboratory equipment, including a large purple balloon, a yellow balloon, and a digital scale. The text "ENZYME SUBSTRATES..." is overlaid in green on the image.


# ENZYME SUBSTRATES...

# ENZYME SUBSTRATES

## CHROMOGENIC SUBSTRATES

Chromogenic, or dye-labelled polysaccharides are useful for the specific measurement of polysaccharide *endo*-hydrolases in crude plant extracts or industrial enzyme preparations. Traditionally, such enzymes have been measured using the native polysaccharide, followed by quantification of the increase in reducing sugar or decrease in viscosity on hydrolysis. Since a range of enzymes, including *endo*- and *exo*-polysaccharidases and glycosidases, act on any given polysaccharide, reducing-sugar methods are not specific. Viscosity reduction methods are specific for *endo*-hydrolase activity, but are tedious to perform, and require specialist equipment. Chromogenic polysaccharide substrates offer the advantages of being specific and sensitive, and can form the basis of accurate, reliable and robust assay procedures.

In the case of insoluble chromogenic substrates the insoluble substrate (gelatinous particles) is depolymerised and solubilised by the action of the *endo*-hydrolase. The reaction is terminated by adding an alkaline solution to stop enzyme activity and the reaction slurry is filtered or centrifuged. Colour in the filtrate or supernatant is measured in a spectrophotometer and the colour intensity is directly related to enzyme activity. Shown below is a schematic representation of the use of AZCL-amylose (Cat. No. I-AZAMY) to measure  $\alpha$ -amylase activity in a sample.



## INSOLUBLE (CROSSLINKED) CHROMOGENIC SUBSTRATES


Cat. No.	Product	For the Measurement of	Recommended Conditions
I-AZAMY	AZCL-Amylose	$\alpha$ -Amylase (Fungal)	Na acetate, 100 mM, pH 4.4
		$\alpha$ -Amylase (Cereal)	Na maleate, 100 mM, pH 6.0
		$\alpha$ -Amylase (Bacterial)	Bis-Tris, 100 mM, pH 7.0
I-AZDAR	AZCL-Arabinan (debranched)	<i>endo</i> -Arabinanase	Na acetate, 50 mM, pH 4.0
I-AZCEL	AZCL-HE-Cellulose	<i>endo</i> -Cellulase ( <i>Trichoderma</i> )	Na acetate, 25 mM, pH 4.5
I-AZXYG	AZCL-Xyloglucan (tamarind)	<i>endo</i> -Cellulase ( <i>Aspergillus</i> )	Na acetate, 25 mM, pH 4.5
I-ACELL	Azo- $\alpha$ -Cellulose		
I-AAVIC	Azo-Avicel		
I-AZBGL	AZCL-Barley $\beta$ -Glucan	<i>endo</i> -Cellulase ( <i>Trichoderma</i> )	Na acetate, 25 mM, pH 4.5
		Lichenase	Na phosphate, 25 mM, pH 6.5
		Malt $\beta$ -Glucanase	Na acetate, 25 mM, pH 4.5
I-AZPAC	AZCL-Pachyman	<i>endo</i> -1,3- $\beta$ -Glucanase	Na acetate, 50 mM, pH 6.0
I-AZCUR	AZCL-Curdlan		



Cat. No.	Product	For the Measurement of	Recommended Conditions
I-AZCHAN	AZCL-Chitosan	<i>endo</i> -Chitosanase	Na acetate, 50 mM, pH 5.0
I-AZDEX	AZCL-Dextran (No. B-512)	<i>endo</i> -1,6- $\alpha$ -Dextranase	Na acetate, 50 mM, pH 5.0
I-AZGLP	AZCL-Galactan (potato)	<i>endo</i> -1,4- $\beta$ -Galactanase	Na acetate, 25 mM, pH 4.3
I-AZGMA	AZCL-Galactomannan (carob)	<i>endo</i> -1,4- $\beta$ -Mannanase	Na acetate, 50 mM, pH 4.5
I-AZCAS	AZCL-Casein	<i>endo</i> -Protease	Na phosphate, 100 mM, pH 7.0
I-AZCOL	AZCL-Collagen		
I-AZPUL	AZCL-Pullulan	Microbial pullulanase	Na acetate, 100 mM, pH 5.0
		Malt limit-dextrinase	Na maleate, 100 mM, pH 5.5
I-AZRHI	AZCL-Rhamnogalacturonan I	Rhamnogalacturonanase	Na acetate, 50 mM, pH 4.5 (or 8)
I-AZXBW	AZCL-Xylan (birchwood)	<i>endo</i> -Xylanase	Na acetate, 25 mM, pH 4.7
I-AZXBE	AZCL-Xylan (beechwood) <sup>NEW</sup>		
I-AZWAX	AZCL-Arabinoxylan (wheat)		

## SOLUBLE CHROMOGENIC SUBSTRATES

In the case of soluble chromogenic substrates, an enzyme sample is incubated with the soluble substrate. The reaction is terminated and high molecular weight, partially hydrolysed fragments are precipitated from solution with an organic solvent such as ethanol or methoxyethanol. Lower molecular weight fragments (products of enzymatic hydrolytic activity) remain in solution. The suspension is mixed thoroughly, centrifuged, and the colour in the supernatant solution is measured in a spectrophotometer. With the aid of a standard curve, the enzyme activity can be determined. Shown below is a schematic representation of the use of Red Starch (Cat. No. S-RSTAR) to measure  $\alpha$ -amylase activity in a sample.



## SOLUBLE CHROMOGENIC SUBSTRATES

Cat. No.	Product	For the Measurement of	Recommended Conditions	Precipitant Solution
S-RDAR	Red Debranched Arabinan (sugar beet)	<i>endo</i> -Arabinanase	Na acetate, 200 mM, pH 4.5	95% v/v Ethanol
S-RSTAR	Red Starch	$\alpha$ -Amylase (Fungal) $\alpha$ -Amylase (Cereal) $\alpha$ -Amylase (Bacterial)	Na malate, 100 mM, pH 5.4 Na malate, 100 mM, pH 5.4 Na maleate, 100 mM, pH 6.5	95% v/v Ethanol
S-ACMCL	Azo-CM-Cellulose (liquid)	<i>endo</i> -Cellulase	Na acetate, 100 mM, pH 4.5	80% v/v Ethanol/ Na Acetate, Zn Acetate, pH 5
S-ACMC	Azo-CM-Cellulose (powder)			
S-AZXG	Azo-Xyloglucan (tamarind)	<i>endo</i> -Cellulase	Na acetate, 100 mM, pH 4.5	95% v/v Ethanol
S-ABG100	Azo-Barley Glucan	<i>endo</i> -Cellulase Lichenase Malt $\beta$ -Glucanase	Na acetate and Na phosphate 40 mM, pH 4.6	80% v/v Methoxyethanol/ Na Acetate, Zn Acetate, pH 5
S-AZFR5	Azo-Fructan	<i>endo</i> -Fructanase	Na acetate, 100 mM, pH 4.5	86% Ethanol/ KOH, 200 mM
S-AZFRXOI	Azo-Fructan plus <i>exo</i> -Inulinase			
S-RPUL	Red Pullulan	Microbial pullulanase Malt limit-dextrinase	Na acetate, 200 mM, pH 5.0 Na acetate, 200 mM, pH 5.0	95% v/v Ethanol
S-ACGLM	Azo-Carob Galactomannan	<i>endo</i> -1,4- $\beta$ -Mannanase	Na acetate, 200 mM, pH 4.0	95% v/v Ethanol
S-AGALP	Azo-Galactan (potato)			
S-AZCAS	Azo-Casein (Sulphanilamide Dyed)	<i>endo</i> -Protease	Na phosphate, 100 mM, pH 7.0	5% Trichloroacetic acid
S-AZRH	AZ-Rhamnogalacturonan	Rhamnogalacturonanase	Na acetate, 100 mM, pH 4.5	95% v/v Ethanol
S-AWAXL	Azo-Wheat Arabinoxylan (liquid)	<i>endo</i> -Xylanase	Na acetate, 50 mM, pH 4.5	95% v/v Ethanol
S-AWAXP	Azo-Wheat Arabinoxylan (powder)			
S-AXBL	Azo-Xylan (birchwood) (liquid)			
S-AXBP	Azo-Xylan (birchwood) (powder)			

## ENZYME TABLET TESTS

Megazyme supplies a range of enzyme tablet tests for ultimate end-user convenience. These products contain the insoluble chromogenic substrates discussed above and from a procedural perspective, operate in the same way. The advantage of enzyme tablet tests is that the need to accurately weigh substrate quantities (and the error associated with this parameter) is removed.


Most enzyme tablet test products are available in pack sizes of 200 and 1000 tablets. They are listed below, separated in terms of the enzyme activities that can be assayed along with the suggested buffer conditions for each test.







Cat. No.	Product	For the Measurement of	Recommended Conditions
T-AMZ	Amylzyme	$\alpha$ -Amylase (Fungal)	Na acetate, 100 mM, pH 4.4
T-AMZBG	Amylzyme BG	$\alpha$ -Amylase (Cereal)	Na maleate, 100 mM, pH 6.0
T-AMZRD	Amylzyme Red	$\alpha$ -Amylase (Bacterial)	Bis-Tris, 100 mM, pH 7.0
T-AMZHY	Amylzyme HY <sup>NEW</sup>		
T-ARZ	Arabinzyme	Arabinase	Na acetate, 50 mM, pH 4.0
T-CAF	Cellzyme AF <sup>NEW</sup>		
T-CCZ	Cellzyme C	endo-Cellulase	Na acetate, 25 mM, pH 4.5
T-CTZ	Cellzyme T		
T-BGZ	Beta-Gluczyme	endo-Cellulase	Na acetate, 25 mM, pH 4.5
		Lichenase	Na phosphate, 25 mM, pH 6.5
		Malt $\beta$ -Glucanase	Na acetate, 25 mM, pH 4.5
T-PAZ	I,3-Beta-Gluczyme	endo-I,3- $\beta$ -Glucanase	Na acetate, 50 mM, pH 6.0
T-CUR	I,3-Beta-Gluczyme HS <sup>NEW</sup>		
T-CHZ	Chitozyme	Chitosanase	Na acetate, 50 mM, pH 5.0
T-DEXT	Alpha-Dextrzyme	endo-I,6- $\alpha$ -Dextranase	Na acetate, 50 mM, pH 5.0
T-LDZ	Limit-Dextrzyme	Microbial pullulanase	Na acetate, 100 mM, pH 5.0
		Malt limit-dextrinase	Na maleate, 100 mM, pH 5.5
T-GLZ	Galactzyme	endo-I,4- $\beta$ -Galactanase	Na acetate, 25 mM, pH 4.3
T-MNZ	Mannzyme	endo-I,4- $\beta$ -Mannanase	Na acetate, 50 mM, pH 4.5
T-PRAK	Protzyme AK	endo-Protease	Na phosphate, 100 mM, pH 7.0
T-PROL	Protzyme OL		
T-RHAM	Rhamnozyme	Rhamnogalacturonanase	Na acetate, 50 mM, pH 4.5 (or 8)
T-XYZ	Xylzyme	endo-Xylanase	Na acetate, 25 mM, pH 4.7
T-XAF	Xylzyme AF <sup>NEW</sup>		
T-XAX	Xylzyme AX		








## COLOURIMETRIC OLIGOSACCHARIDES

The colourimetric substrates supplied by Megazyme are based on defined oligosaccharides which are covalently linked to a colourimetric or fluorimetric moiety through the reducing D-glucosyl residue of the oligosaccharide. In the measurement of  $\alpha$ -amylase using Ceralpha reagent, the substrate is composed of end-blocked 4-nitrophenyl- $\alpha$ -maltoheptaoside in the presence of an excess quantity of thermostable  $\alpha$ -glucosidase. When the substrate is cleaved by  $\alpha$ -amylase, the  $\alpha$ -glucosidase removes the remaining D-glucosyl residues releasing free 4-nitrophenol, which in the presence of an alkaline solution is converted to the yellow phenolate ion. The principle is depicted below.



Colourimetric oligosaccharides that do not contain a blocking group (e.g. 2-chloro-4-nitrophenol-cello-oligosaccharides) cannot be used to measure *endo*-acting hydrolytic in the presence of *exo*-acting enzymes but can be employed for the assay of pure enzymes. These substrates also find use in enzyme research employing HPLC analysis where their chromogenic moiety allows for quantification by UV detection.

Cat. No.	Product	
<b>2-Chloro-4-Nitrophenyl-Cello-Oligosaccharides</b>		
O-CPNPG2	DP2	2-Chloro-4-nitrophenyl- $\beta$ -cellobioside <sup>NEW</sup>
O-CPNPG3	DP3	2-Chloro-4-nitrophenyl- $\beta$ -cellotrioside <sup>NEW</sup>
O-CPNPG4	DP4	2-Chloro-4-nitrophenyl- $\beta$ -cellotetraoside <sup>NEW</sup>
O-CPNPG5	DP5	2-Chloro-4-nitrophenyl- $\beta$ -cellopentaoside <sup>NEW</sup>
		
<b>4-Methylumbelliferyl-beta-Cello-Oligosaccharides</b>		
O-4MUG2	DP2	4-Methylumbelliferyl- $\beta$ -cellobioside <sup>NEW</sup>
O-4MUG3	DP3	4-Methylumbelliferyl- $\beta$ -cellotrioside <sup>NEW</sup>
O-4MUG4	DP4	4-Methylumbelliferyl- $\beta$ -cellotetraoside <sup>NEW</sup>
O-4MUG5	DP5	4-Methylumbelliferyl- $\beta$ -cellopentaoside <sup>NEW</sup>
		
O-B4MUG3	Blocked 4-methylumbelliferyl- $\beta$ -cellotrioside <sup>NEW</sup>	
		
O-BNAPG3	Blocked $\beta$ -naphthyl- $\beta$ -cellotrioside <sup>NEW</sup>	
		

Cat. No.	Product		
<b>4-Nitrophenyl-β-Xylo-Oligosaccharides</b>			
O-PNPX2	DP2	4-Nitrophenyl-β-xylobioside <sup>NEW</sup>	
O-PNPX3	DP3	4-Nitrophenyl-β-xylotrioside <sup>NEW</sup>	
O-ONPX2	DP2	2-Nitrophenyl-β-xylobioside <sup>NEW</sup>	
<b>4-Methylumbelliferyl-β-Xylo-Oligosaccharides</b>			
O-4MUX2	DP2	4-Methylumbelliferyl-β-xylobioside <sup>NEW</sup>	
O-4MUX3	DP3	4-Methylumbelliferyl-β-xylotrioside <sup>NEW</sup>	
<b>2-Chloro-4-Nitrophenyl-α-Manno Oligosaccharides</b>			
O-CPNPAM2	DP2	2-Chloro-4-nitrophenyl-α-mannobioside <sup>NEW</sup>	
O-CPNPAM3	DP3	2-Chloro-4-nitrophenyl-α-mannotrioside <sup>NEW</sup>	
O-CPNPAM4	DP4	2-Chloro-4-nitrophenyl-α-mannotetraoside <sup>NEW</sup>	
<b>4-Methylumbelliferyl-α-Manno-Oligosaccharides</b>			
O-4MUAM2	DP2	4-Methylumbelliferyl-α-mannobioside <sup>NEW</sup>	
O-4MUAM3	DP3	4-Methylumbelliferyl-α-mannotrioside <sup>NEW</sup>	
O-4MUAM4	DP4	4-Methylumbelliferyl-α-mannotetraoside <sup>NEW</sup>	
<b>Nitrophenyl-Malto-Oligosaccharides</b>			
O-BPNPC7	Blocked 4-nitrophenyl-α-maltoheptaoside		
O-PNPC3	4-Nitrophenyl-β-maltotrioside <sup>NEW</sup>		
O-NAPC3	β-Naphthyl-β-maltotrioside <sup>NEW</sup>		
O-PNPL	4-Nitrophenyl β-lactoside <sup>NEW</sup>		





CARBOHYDRATES...






# CARBOHYDRATES

Megazyme provides high purity oligosaccharides and polysaccharides for research purposes. Although in many cases the measurement of enzymatic activity is greatly simplified through the use of chromogenic / colourimetric substrates, activity on the native substrates can also be investigated. Megazyme's oligosaccharides and polysaccharides are also used all over the world for investigations into the mechanism of action of enzymes. The common oligosaccharide families (cello-, xyl-, manno- etc.) are generally available in DP2-6 and are accompanied by a data sheet describing purity by HPLC/TLC as appropriate. The most common polysaccharides are available in a range of viscosity grades that vary in average molecular weight, degree of branching on the polysaccharide backbone, or both. Accompanying data sheets outline relevant physicochemical data.




## OLIGOSACCHARIDES


Cat. No.	Product	
Acetyl-Chito-Oligosaccharides		
		
O-CHI2	DP2	Diacetyl-chitobiose
O-CHI3	DP3	Triacetyl-chitotriose
O-CHI4	DP4	Tetraacetyl-chitotetraose
O-CHI5	DP5	Pentaacetyl-chitopentaose
O-CHI6	DP6	Hexaacetyl-chitohexaose

### Aldouronic Acids (from xylan)

O-AMX	Aldouronic acid mixture	
O-AMXR	Aldouronic acid mixture (borohydride reduced)	

### Amylosaccharides (mixed-linkage)


	
O-IPAN	Isopanose <sup>NEW</sup>
	
O-GMT	6'-alpha-D-Glucosyl-maltotriose
	
O-GMH	6'-alpha-D-Glucosyl-maltotriosyl-maltotriose

Cat. No.	Product
	
O-MTMT	6'- $\alpha$ -D-maltotriosyl-maltotriose
O-MTMTRD	6'- $\alpha$ -D-maltotriosyl-maltotriose (BH <sub>4</sub> reduced) <sup>NEW</sup>


### 1,5-α-L-Arabetriose-Oligosaccharides

O-ABI	DP2	Arabinobiose (syrup)
O-ATR	DP3	Arabinotriose (syrup)
O-ATE	DP4	Arabinotetraose (syrup)
O-APE	DP5	Arabinopentaose (syrup)
O-AHE	DP6	Arabinohexaose (powder)
O-AHP	DP7	Arabinoheptaose (powder)
O-AOC	DP8	Arabino-octaose (powder)

### Cello-Oligosaccharides


		
O-CTR	DP3	Cellotriose
O-CTE	DP4	Cellotetraose
O-CPE	DP5	Cellopentaose
O-CHE	DP6	Cellohexaose

### Borohydride Reduced Cello-Oligosaccharides


		
O-CTRRD	DP3	1,4-β-D-Cellotriitol <sup>NEW</sup>
O-CTERD	DP4	1,4-β-D-Cellotetraitol <sup>NEW</sup>
O-CPERD	DP5	1,4-β-D-Cellopentaitol <sup>NEW</sup>
O-CHERD	DP6	1,4-β-D-Cellohexaitol <sup>NEW</sup>

**Cat. No. Product**


**1,3:1,4 β-Gluco-Oligosaccharides**




O-BGTRIA 1,3:1,4-β-Glucotriose A




O-BGTRIB 1,3:1,4-β-Glucotriose B



O-BGTETA 1,3:1,4-β-Glucotetraose A




O-BGTETB 1,3:1,4-β-Glucotetraose B




O-BGTETC 1,3:1,4-β-Glucotetraose C

**Fructo-Oligosaccharides**




O-KTR DP3 1-Kestose  
O-KTE DP4 1,1-Kestotetraose  
O-KPE DP5 1,1,1-Kestopentaose

**Galacto-Manno-Oligosaccharides**




O-GM3 6'-α-D-Galactosyl-mannotriose

**Cat. No. Product**




O-GMM3 6'-α-D-Galactosyl-mannobiose + mannotriose




O-GGM5 6³,6⁴-α-D-Galactosyl-mannopentaose

**Galacto-Oligosaccharides**




O-GBI 1,3(4)-β-D-Galactobiose (purity > 95%)

**Galactosyl-Sucrose Oligosaccharides**




O-VER Verbascose

**Gluco-Manno-Oligosaccharides**




O-GMMBI 1,4-β-D-Glucosyl-D-Mannose plus 1,4-β-D-Mannobiose<sup>NEW</sup>




O-GMMTR 1,4-β-D-Glucosyl-D-Mannobiose and 1,4-β-D-Cellobiosyl-D-Mannose<sup>NEW</sup>


## Cat. No. Product

**1,3-β-D-Gluco-Oligosaccharides**


O-LAM2	DP2	Laminaribiose
O-LAM3	DP3	Laminaritriose
O-LAM4	DP4	Laminaritetraose
O-LAM5	DP5	Laminaripentaose
O-LAM6	DP6	Laminarihexaose

**1,4-β-D-Manno-Oligosaccharides**


O-MBI	DP2	Mannobiose
O-MTR	DP3	Mannotriose
O-MTE	DP4	Mannotetraose
O-MPE	DP5	Mannopentaose
O-MHE	DP6	Mannohexaose

**Xyloglucan Derived Oligosaccharides**


O-IPRM	Isoprimeverose (xyloglucan derived)
--------	-------------------------------------



O-XCBIR	Xylosyl-cellobiose (NaBH <sub>4</sub> reduced)
---------	--




O-X3G4	Xyloglucan heptasaccharide (X <sub>3</sub> Glc <sub>4</sub> )
--------	---




O-X3G4R	Heptasaccharide (X <sub>3</sub> Glc <sub>4</sub> ; NaBH <sub>4</sub> reduced)
O-XGHON	Xyloglucan (Hepta + octa + nona-saccharides)
O-XGHDP	Higher DP xyloglucan oligosaccharides

## Cat. No. Product


**1,4-β-D-Xylo-Oligosaccharides**

O-XBI	DP2	Xylobiose
O-XTR	DP3	Xylotriose
O-XTE	DP4	Xylotetraose
O-XPE	DP5	Xylopentaose
O-XHE	DP6	Xylohexaose


O-XBIRD	DP2	Xylobiose (borohydride reduced) <sup>NEW</sup>
O-XTRRD	DP3	Xylotriose (borohydride reduced) <sup>NEW</sup>

**Arabino-Xylo-Oligosaccharides**



O-AX2	3 <sup>2</sup> -α-L-Arabinofuranosyl-xylobiose <sup>NEW</sup>
-------	---



O-AX3	2 <sup>3</sup> -α-L-Arabinofuranosyl-xylotriose <sup>NEW</sup>
-------	--



O-A2X3	2 <sup>3</sup> ,3 <sup>3</sup> -di-α-L-Arabinofuranosyl-xylotriose <sup>NEW</sup>
--------	---



O-AX3MIX	2 <sup>3</sup> -α-L-Arabinofuranosyl-xylotriose <sup>NEW</sup> plus 3 <sup>3</sup> -α-L-Arabinofuranosyl-xylotriose
----------	--

## POLYSACCHARIDES

Cat. No.	Product
P-ARAB	Arabinan (sugar beet)
P-DBAR	Debranched Arabinan (sugar beet)
P-LARB	Linear 1,5- $\alpha$ -L-Arabinan (sugar beet)
P-CMLA	CM-Linear 1,5- $\alpha$ -L-Arabinan (sugar beet)
P-ARGAL	Arabinogalactan (larch wood)
P-RAXY	Arabinoxylan (rye flour)
P-WAXYL	Arabinoxylan (wheat flour; low viscosity)
P-WAXYM	Arabinoxylan (wheat flour; medium viscosity)
P-WAXYH	Arabinoxylan (wheat flour; high viscosity)
P-WAXYI	Arabinoxylan (wheat flour; insoluble)
P-ADWAX22	Arabinoxylan (acid debranched; 22% ara) <sup>NEW</sup>
P-ADWAX26	Arabinoxylan (acid debranched; 26% ara) <sup>NEW</sup>
P-EDWAX30	Arabinoxylan (enzyme debranched; 30% ara) <sup>NEW</sup>
P-BGBL	Beta-Glucan (barley; low viscosity)
P-BGBM	Beta-Glucan (barley; medium viscosity)
P-BGBH	Beta-Glucan (barley; high viscosity)
P-BGOM	Beta-Glucan (oat; medium viscosity)
P-BGOH	Beta-Glucan (oat; high viscosity)
P-BGYST	Beta Glucan (yeast; alkali soluble)
P-BGCFA	Beta-Glucan CFA standard
P-MWBGs	Beta-Glucan MW standard
P-BLDX	Beta-Limit Dextrin (10 g)
P-BLDX50	Beta-Limit Dextrin (50 g)
P-CMC4M	Carboxymethyl Cellulose 4M
P-CHITIN	Chitin (colloidal) <sup>NEW</sup>
P-CHITOSAN	Chitosan <sup>NEW</sup>
P-CURDL	Curdlan
P-CMCUR	CM-Curdlan

P-FOS28	Fructooligosaccharides (DP2-8) <sup>NEW</sup>
P-GALLU	Galactan (lupin)
P-GALPOT	Galactan (potato)
P-GALML	Galactomannan (carob; low viscosity)
P-GALMH	Galactomannan (carob; high viscosity)
P-GGMMV	Galactomannan (guar; medium viscosity)
P-GGMMHV	Galactomannan (guar; high viscosity)
P-GGM2I	Galactomannan (guar; galactose depleted; 21% gal)
P-GGM28	Galactomannan (guar; galactose depleted; 28% gal)
P-GLCML	Glucomannan (konjac; low viscosity)
P-GLCMH	Glucomannan (konjac; high viscosity)
P-INUL	Inulin (DP2-60) <sup>NEW</sup>
P-LICHN	Lichenan (icelandic moss)
P-MANIV	Mannan (ivory nut)
P-MANCB	Mannan (1,4- $\beta$ -D-mannan)
P-PACHY	Pachyman (1,3- $\beta$ -D-glucan)
P-CMPAC	CM-Pachyman
P-PGALU	Pectic Galactan (lupin)
P-PGAPT	Pectic Galactan (potato)
P-PGACT	Polygalacturonic Acid (PGA)
P-PULLN	Pullulan
P-PULLBH	Pullulan (NaBH <sub>4</sub> reduced)
P-RHAMI	Rhamnogalacturonan I (potato)
P-RHAGN	Rhamnogalacturonan (soy bean)
P-SCLER	Scleroglucan <sup>NEW</sup>
P-XYGLN	Xyloglucan (tamarind)
P-XYLNBE	Xylan (Beechwood; purified) <sup>NEW</sup>







# EQUIPMENT...

## EQUIPMENT

Megazyme provides a number of equipment items for use with its diagnostic applications. The wine industry in particular is well catered for with a range of spectrophotometers, from the Megaquant™ Meter which is suitable for any home brewer right up to the fully automated Chemwell® 2910 which is sufficient for a large winery or wine analysis laboratory.

### Cat. No. D-CHEM2910

- ◆ Suitable for large wineries / laboratories
- ◆ Automated assay formats (200 tests per hour)
- ◆ Pre-programmed with Megazyme test settings
- ◆ Full range of tests available
- ◆ Simple assay procedures


### ChemWell®2910



### Cat. No. D-CHEMT

- ◆ Suitable for medium wineries / laboratories
- ◆ Automated assay formats (100 tests per hour)
- ◆ Pre-programmed with Megazyme test settings
- ◆ Full range of tests available
- ◆ Simple assay procedures

### ChemWell®-T



### Cat. No. D-STATFAX

- ◆ Suitable for small wineries / laboratories
- ◆ Manual assay formats
- ◆ Pre-programmed with Megazyme test settings
- ◆ Full range of tests available
- ◆ Simplified data analysis
- ◆ Reagent stability > 2 years

### Stat Fax® 4500



### Cat. No. D-FRGLMQ or D-LMALMQ

- ◆ Small and portable
- ◆ Simple to use
- ◆ Compatible with K-FRGLMQ and K-LMALMQ
- ◆ Extremely affordable

### Megaquant™





Cat. No.	Product
<b>EQUIPMENT</b>	
D-CHEM2910	ChemWell®2910 Automated EIA & Chemistry Analyser <sup>NEW</sup>
D-CHEMT	ChemWell®-T Automated Chemistry Analyser <sup>NEW</sup>
D-STATFAX	Stat Fax® 4500 Chemistry Analyser (Spectrophotometer) <sup>NEW</sup>
D-SFTUBE	Tubes for Stat Fax® 4500 Chemistry Analyser <sup>NEW</sup>
D-FRGLMQ	MegaQuant™ Meter plus D-Fructose & D-Glucose Reagents
D-LMALMQ	MegaQuant™ Meter plus L-Malic Acid Reagents
D-MQTUB	Tubes for MegaQuant™ Meter (24 tubes)
D-IBMKIII	Megazyme Incubation Bath MK III
D-INTDFB	Water Bath for Integrated Total Dietary Fibre Procedure <sup>NEW</sup>
<b>BOOK</b>	
D-ADFTB	Advanced Dietary Fibre Technology Book
<b>GENERAL CHEMICALS</b>	
G-AMBOH	Amberlite FPA OH <sup>-</sup> Ion Exchange Resin
G-AMBH	Ambersep 200 H <sup>+</sup> Ion Exchange Resin
G-CELITE	Celite
G-LCYST200	L-Cysteine Hydrochloride Monohydrate

Cat. No.	Product
<b>COFACTORS &amp; STAINS</b>	
C-ATP	Adenosine 5'-triphosphate
C-CLFR	Calcofluor fluorescent stain
C-COA500	Coenzyme A (trilithium salt)
C-NAD	β-Nicotinamide adenine di-nucleotide
C-NADP	β-Nicotinamide adenine di-nucleotide phosphate
C-NADH	β-Nicotinamide adenine di-nucleotide reduced salt
<b>BUFFERS</b>	
B-BISTRIS250	BIS-TRIS Buffer Salt
B-CAPS200	CAPS Buffer Salt
B-CAPSO250	CAPSO Buffer Salt
B-GLYGLY250	Glycylglycine Buffer Salt
B-HEPES250	HEPES Buffer Salt
B-MES250	MES Monohydrate Buffer Salt
B-MOPS250	MOPS Buffer Salt
B-PIPES250	PIPES Buffer Salt
B-TRIS500	TRIS Buffer Salt
<b>LECTINS</b>	
L-CONA	Concanavalin A





The screenshot shows the Megazyme website interface. At the top, there is a navigation bar with links: Quick Order, Contact Us, Worldwide Distributors, Cart, EUR, and Login/Register. Below this is a green header with the Megazyme logo on the left, the tagline 'Setting New Standards in Test Technology' in the center, and a search bar on the right. A secondary navigation bar contains links: Home, Products, Services, Resources, Media Centre, About Us, and Help.

The main content area features a large banner image of a laboratory with the word 'Innovation' overlaid in large green letters. Below the banner, there are three main sections:

- Select an Industry:** Includes icons and links for Fermentation, Food, Brewing, and Feed.
- Latest News:** Features a news item titled 'Megazyme is named overall winner at the SFA National Small Business Awards 2013' dated Mar 08, 2013, with a 'read more' link.
- New Products:** Lists products such as 'endo-Cellulase Assay Kit' and 'L-Arabinose / D-Galactose Assay Kit', with a 'View All' link.

At the bottom of the website, there is a footer with links for Legal, SSL/Privacy, Cookies, Terms & Conditions, and Site Map. It also includes contact information for Megazyme International Ireland and a Newsletter Sign Up form with fields for email, first name, and last name, and a 'sign up' button.

[www.megazyme.com](http://www.megazyme.com)

**Megazyme also has a team of International Distributors in over 45 countries. Contact details can be found online at [www.megazyme.com](http://www.megazyme.com)**



**Megazyme,  
Bray Business Park, Bray,  
Co. Wicklow, A98 YV29, Ireland.**

**T + 353 1 286 1220**

**Megazyme Inc.,  
55 East Monroe Street, Suite 3800  
21st Floor, Chicago, Illinois, 60603, USA.**

**T 312 212 4361**

**E [cs@megazyme.com](mailto:cs@megazyme.com)**

**ISO 9001:2008 Registered**

**[www.megazyme.com](http://www.megazyme.com)**



**Cover image hand-painted by Lisa  
McCleary.  
[lmccleary2@gmail.com](mailto:lmccleary2@gmail.com)  
[mcc-art.tumblr.com](http://mcc-art.tumblr.com)**

*Lisa McCleary*





Lisa McCreary



